

STRETCHER REPLACEMENT MATTRESSES

Mattress Compatible with Stryker XPS

DSC-9577F-0518Q74



- Designed to increase patient comfort while reducing transfer gap
- Heavy-duty vinyl is highly resistant to tearing
- Compatible with Stryker Power-PRO XT Model 6500/6506, Stryker Power-PRO TL Model 6550, and Stryker Performance-PRO XT Model 6085/6086
- Must be used in conjunction with Stryker XPS side rails

Universal Stretcher Replacement Mattress

DSC-9577-0518H74



- The 'bolster' styling cradles the patient for enhanced comfort
- Compatible with all non-XPS Stryker EMS stretchers
- Heavy-duty vinyl is highly resistant to tearing

DAVID SCOTT
COMPANY