




Available in various thicknesses and foam configurations, our pads ensure the perfect fit for your needs. Our replacement surgical table pads have a lifespan of 3 to 5 years, offering reliable performance in the operating room. For any questions or to inquire about custom pads, please contact our sales team.



SHEERCARE SERIES (LEVEL I)

This pad is constructed of high density base foam surrounded by a Fusion III® ultra soft 4-way stretch covering. Our SheerGuard Series offers an economical alternative without compromising the comfort and safety of your patients.

[High Density Base Foam](#) → [Fusion III® Ultra Soft Covering](#)



ULTRACARE SERIES (LEVEL II)

These pad sets are constructed with high density base memory foam, then topped with a ½” layer of Visco memory foam and surrounded by a Fusion III® ultra soft 4-way stretch covering.

[High Density Base Foam](#) → [Visco Memory Foam](#) → [Fusion III® Ultra Soft Covering](#)



COMBO SERIES (LEVEL III)

Our unique “tri-level” design is constructed of a high-density supportive base foam then topped with a ½” layer of doughy soft cradling Visco memory foam. A layer of Blue Diamond® Gel is placed just below the Fusion III® ultra soft 4-way stretch covering to ensure your patients the maximum protection available for pressure management.

[High Density Base Foam](#) → [Visco Memory Foam](#) → [Blue Diamond Gel®](#) → [Fusion III® Ultra Soft Covering](#)