

# ProPette™ Digital

## *Electronic Pipette Controller*



- ⦿ Enhanced flow control - Max aspirate and dispense rates
- ⦿ Independent aspiration and dispense speed settings
- ⦿ Pressure sensitive triggers for even more control
- ⦿ Improved QuickStand™ - Locks in place with a click
- ⦿ USB-C charge port - No searching for specific chargers



ProPette™ Digital Pipette Controllers add even more convenience and functionality to the popular ProPette™ Line.

New buttons and an enhanced digital display allow users to set their max aspirate and dispense rates on the fly. Paired with ProPette™ pressure-sensitive triggers, users can enjoy the most precise flow controls possible. The display also keeps constant track of the battery level, which can be charged via the convenient new USB-C port.

# ProPette™ Digital Electronic Pipette Controller



**Digital Display**  
Battery Level  
& Max Flow Rates



**Flow Controls**  
Pressure-Sensitive  
Control Triggers



**USB-C Charging**  
Universal Charge Port  
w/ Included Cable

ProPette™ Controllers from MTC Bio™ are lightweight (only 6.3oz/170g) and designed to fit comfortably in either the right or left hand. The natural feel of the ergonomic body design, the computer-balanced weight distribution, and the low-pressure fingertip control button all combine to provide a more comfortable pipetting experience. The ProPette™ controller includes the proprietary Quickstand™ retractable stand, which allows the instrument to be set down on the bench without risk of contamination.



Cat No.	Description	Qty
P6290	ProPette™ Digital Controller, Digital Display and Controls	1 / Each
P6080-BA	Replacement Battery, Lithium, 3.7V, 750mAh for Pipette Controller	1 / Each
P6080-FI	Replacement Filters, 0.8µm for Pipette Controller	5 / Pack
P6080-NA	Silicone Nosepiece Adapter for Pipette Controller	1 / Each
P6080-PF	LiquiLock™ Premium Pipette Filter, Protects from Liquid Overflow	5 / Pack