



# BIGGER, BETTER BenchBin 3L™

- Twice the capacity of the original BenchBin™
  - Maintains a compact footprint on the bench
- Safely dispose of tips, tubes, gloves, etc.
- Compliant with OSHA 29 CFR 1910.1030
- Secures and covers hazardous waste

The new BenchBin 3L™ is a larger capacity biohazard bin boasting all of the key features that made the original BenchBin™ an instant favorite for benchtop waste disposal.

To close the lid, simply push it down and feel it click into place. A second push opens it up again.

BenchBin 3L™ has a stable base that never risks tipping and delivers an overall better user experience than traditional wire rack stands.

BenchBin 3L™ prioritizes a compact footprint (6×5 inch) yet contains twice the capacity of the original BenchBin™

BenchBin 3L™ - A8403B



3L Capacity  
6×5" Footprint



BenchBin™ - A8401B



1.5L Capacity  
7×5" Footprint





# BIGGER, BETTER BenchBin 3L™



Cat. No.	Description	Qty
<b>BenchBin 3L™ and Autoclave Bags</b>		
A8403B	BenchBin 3L™ Benchtop Biohazard Bin, 3L Capacity, Spring-Loaded Lid, 6×5×8.5in. <i>Fits 10×12.5in. Autoclave Bags (A9006R)</i>	1
A9006R	Autoclave Bags, 10×12.5in. Biohazard, Printed, Marking Area, <i>Fits BenchBin 3L™ (A8403B)</i>	100
A8403-SK	BenchBin 3L™ Starter Kit: 1 BenchBin 3L™ and 400 Red 10×12.5in. Autoclave Bags	1
<b>Original BenchBin™ and Autoclave Bags</b>		
A8401B	BenchBin™ Benchtop Biohazard Bin, 1.5L Capacity, Spring-Loaded Lid, 7×5×7in. <i>Fits 8.5×11in. Autoclave Bags (A9004R)</i>	1
A9004R	Autoclave Bags, 10×12.5in. Biohazard, Printed, Marking Area, <i>Fits Original BenchBin™ (A8401B)</i>	100
A8401-SK	BenchBin™ Starter Kit: 1 BenchBin™ and 400 Red 8.5×11in. Autoclave Bags	1