


# PCR Double Decker™ Hinged-Lid Freezer Boxes

- Fits 0.2mL / 0.1mL PCR consumables
  - Tubes, strips, or plates
- Two levels of 14 × 14 (196) wells
- Max capacity: 392 × PCR Tubes
- Standard dimensions (5¼ × 2 inch)
- Hinged lid stays attached
- Drain slots on bottom



## Product Specifications

**Material:** Coated Cardboard

**Dimensions:** 5¼ × 5¼ × 2 inch


133 × 133 × 51mm

## Content Specifications

**Compatible Tubes:** PCR, 0.1 / 0.2mL

**Tube Capacity:** 392×0.2mL (2×196)

**Strip Capacity:** 40×8-strips, 28×12-strips



# PCR Double Decker™ Hinged-Lid Freezer Boxes


This unique new freezer box has been designed to **maximize storage capacity for PCR tubes.**

With 196 wells on each level, it can accommodate 392 individual 0.1mL or 0.2mL tubes. For strip tubes, the box handles 40 strips of 8 or 28 strips of 12.

The top level is a cardboard tray with two finger grips that allow easy access to the lower level. The two level system doubles the capacity while keeping the total box height under 2 inches (51mm) and assuring a good fit in standard height stainless steel racks.



Box and Lower Level



Removable Upper Tray

CAT. NO.	DESCRIPTION	Quantity
R2392	Double Decker™ Cardboard PCR Box, 5¼ × 5¼ × 2 inch, 392 × PCR Tubes	Pack of 5