



Buried in STYROFOAM RACKS?

Introducing a *Green Alternative!*
Eco-Plastic Racks for 15 and 50mL Tubes




Same great tubes in a
new eco-friendly package

- Completely recyclable fold-flat racks
- Reduce waste volume by 65%
- 25 tubes per rack sealed in a sterile bag
- More durable and reusable than foam



FOLDS FLAT!
65% Less Waste Volume



Introducing a **Green Alternative!**


Eco-Plastic Racks for 15 and 50mL Tubes

In the same way that styrofoam food and drink containers are being phased out and fully banned in some locations, research laboratories are also working toward reducing their environmental footprint and increasing sustainability.

Polystyrene foam tube racks have an overwhelming environmental impact and can **take several hundred years to decompose***



**Recyclable,
Reusable, and
Autoclavable**



Our new Eco-Racks contain the same high quality centrifuge tubes that you've come to expect from MTC Bio, but with a new focus on sustainability. With an easy-turn cap that goes on smoothly and seals completely, our tubes are guaranteed to remain leakproof to 95kPa/14PSI.

Made of high clarity polypropylene with printed graduations and marking spots, these tubes are made to standard dimensions to fit all popular equipment.

They are supplied sterile and are RNase/DNase/Pyrogen free.

Cat No.	Description	Qty
C2607	15mL PP Tube (17×121mm) Flat Screw Cap, Eco-Plastic Racks of 25 in Sterile Bags	500
C2608	50mL PP Tube (29×117mm) Flat Screw Cap, Eco-Plastic Racks of 25 in Sterile Bags	500