

## GTR-2-AS / Adjustable Speed Tube Rocker



*GTR-2-AS adjustable speed two-tier tube rocker*

*When one platform isn't sufficient, opt for the space-efficient design of Globe Scientific's GTR-2-AS!*

The GTR-2-AS adjustable-speed tube rocker from Globe Scientific is an excellent choice for essential laboratory tube mixing tasks. Its rocking motion, adjustable speed ranging from 0 to 30 RPM, and 30° angle ensure optimal performance. The included silicone inserts accommodate tubes from 0.5mL to 10mL, allowing users to simply place the tubes, switch it on, and leave it to operate.

**Features:**

- Adds an additional platform tier for double the capacity
- Easy plug-and-play functionality
- Compact design requires minimal bench space
- Accommodates tubes from 0.5mL to 10mL (50mL with optional insert)



*If a double platform is more than you need, ask about the Globe Scientific GTR-AS single tier tube rocker.*

## Tube Rocker



Radial dial for simple speed adjustment from 0 to 30rpm



The included ten position silicone inserts are compatible with 5-10mL specimen collection tubes



GTR-AS shown with optional silicone insert for 50mL centrifuge tubes (available separately)

Specifications	GTR-AS
Working Area (each tier)	190.5mm W x 127mm D (7.5in W x 5in D)
Maximum Load	1kg / (2.2lbs)
Maximum Tilt Angle	30°
Voltage	100-240VAC, 50-60Hz
Power	15W
Speed	0-30rpm adjustable speed
Dimensions (WxDxH)	295x160x165mm / (11.6 x 6.3 x 6.5in)
Weight	2kg / (4.4lbs)
Operating Temperature	5-50°C
Permissible Relative Humidity	80%RH (non-condensing)
Warranty	2 Years

### ORDERING INFORMATION

Item #	Description	Unit
GTR-2-AS	Double Tier Tube Rocker, Adjustable Speed, 0-30rpm, 100-240VAC, 50-60Hz	Each
GTR-AS-10X5-10ML	10-Position Silicone Insert for 5-10mL specimen collection tubes, for GTR-AS and GTR-2-AS	Each
GTR-AS-4X50ML	4-Position Silicone Insert, for 50mL centrifuge tubes, for GTR-AS and GTR-2-AS	Each



GTR-AS-10X5-10ML

Replacement inserts for 5-10mL and 50mL tubes are available



GTR-AS-4X50ML

