

## ***Tubes & Vials with Separately Packed Caps*** **Detached Cap Economy Packs**

- 🔗 Facilitates filling and handling
- 🔗 Ideal for working with large quantities
- 🔗 High quality polypropylene tubes
- 🔗 Caps packaged in separate bags



MTC's most popular tubes and vials are now available in new bulk packaging configurations. The tubes and caps are packaged in separate bags, thereby expediting any process that normally includes removing the attached cap prior to filling the tube.



CAT. NO.	DESCRIPTION - 10 Sterile Bags of 100 Tubes	QTY
<b>C1811-SC</b>	Transport Tube, Sterile, 5mL, 16×60mm, with Separate Screw Caps,	1000
<b>C1812-SC</b>	Transport Tube, Sterile, 10mL, 16×80mm, with Separate Screw Caps,	1000

CAT. NO.	DESCRIPTION - 10 Non-Sterile Bags of 100 Tubes	QTY
<b>C1813-NS</b>	SureSeal™ Transport Tube, 5mL, 16×56mm, with Separate O-Ring Caps, Non-sterile	1000
<b>C1814-NS</b>	SureSeal™ Transport Tube, 10mL, 16×92mm, with Separate O-Ring Caps, Non-sterile	1000

CAT. NO.	DESCRIPTION - Non-Sterile Bulk Pack	QTY
<b>C2601-B</b>	Centrifuge Tube, 15mL, 17×121mm, with Separate Caps, Non-sterile	500
<b>C2603-B</b>	Centrifuge Tube, 50mL, 29×117mm, with Separate Caps, Non-sterile	500

# Tubes & Vials with Separately Packed Caps Detached Cap Economy Packs

## Transport Tubes with Screw Cap

*Sterile Bags of 100*



### Specifications

<b>Materials, Tube:</b>	Polypropylene
<b>Materials, Caps:</b>	HD Polyethylene
<b>Temperature:</b>	-86°C - 121°C
<b>Leak-proof:</b>	to 95kPa
<b>Dimensions:</b>	5mL: 16×60mm 10mL: 16×80mm

## SureSeal™ Transport Tubes w/ Silicone O-Ring Screw Cap

*Non-Sterile Bags of 100*




### Specifications

<b>Materials, Tube:</b>	Polypropylene
<b>Materials, Caps:</b>	Polypropylene Silicone O-Ring
<b>Temperature:</b>	-86°C - 121°C
<b>Leak-proof:</b>	to 95kPa
<b>Dimensions:</b>	5mL: 16×56mm 10mL: 16×92mm

## Centrifuge Tubes 15mL and 50mL Tubes

*Non-Sterile Bulk Bags*



### Specifications

<b>Materials, Tube:</b>	Polypropylene
<b>Materials, Caps:</b>	HD Polyethylene
<b>Temperature:</b>	-86°C - 121°C
<b>Leak-proof:</b>	to 95kPa
<b>Max RCF:</b>	15mL: 12,000×g 50mL: 15,000×g
<b>Dimensions:</b>	15mL: 17×121mm 50mL: 29×117mm