




Multichannel Reagent Reservoirs

8 and 12 Channel for Manual or Robotic Pipetting

-  8 or 12 Channel Reservoirs
-  Automation/Multichannel Pipette Compatible
-  Sterile, individually wrapped



These economical disposable reservoirs are made of rigid polystyrene. Ideal for repetitive and robotic liquid transfer applications, they are also designed for use with multichannel pipettes.

8-Channel Reagent Reservoir
P8208 - 8×6.0mL - 48mL Total

12-Channel Reagent Reservoir
P8212 - 12×4.0mL - 48mL Total




Material
Polystyrene

Dimensions
84×127×14mm

Capacity
48mL Total

Sterilization
Electron Beam



Cat No.	Description	Qty
P8208	8-Channel Reagent Reservoir, 48mL (8×6.0mL), E-Beam Sterilized, Individually Wrapped	100 / cs
P8212	12-Channel Reagent Reservoir, 48mL (12×4.0mL), E-Beam Sterilized, Individually Wrapped	100 / cs