

# DIGIDOP "CLASSIC" DD-700

Newman  Medical 



FETAL



VASCULAR



HR DISPLAY



ALKALINE



American-Made  
Quality



Includes labor, parts,  
and accidental  
breakage



Includes cords and  
probes



Waterproof option

## HR Display Doppler System

Don't settle for less than the best with the DigiDop Classic Doppler. This handheld device features Digital Sound Clarity™, ensuring exceptional clarity and accurate readings during examinations. The Classic's standard size provides a comfortable grip and allows for easy portability within your practice. Furthermore, the DigiDop Classic boasts a long battery life, eliminating downtime during your workday.

The DigiDop Classic is built to last. Manufactured in the USA with high-quality materials, each DigiDop undergoes rigorous testing by Newman Medical to guarantee unwavering performance. Invest in the DigiDop Classic and experience the confidence that comes with a trusted and reliable Doppler system.

## DigiDop Doppler - The Clear Choice in Diagnostics

### Unmatched Clarity, Unwavering Support:

- **Hear the Difference:** Digital Sound Clarity DSC™ amplifies the signal up to 6x for exceptional diagnostic accuracy.
- **American-Made Advantage:** Designed, manufactured, and serviced in the U.S. for guaranteed quality.
- **Built to Endure:** Unbreakable probe holder ensures longevity and reliability.
- **Seamless Flexibility:** All probes are interchangeable for maximum efficiency.

### Confidence-Inspiring Warranty:

- **5-Year Worry-Free Coverage:** Our honest warranty covers all parts and labor, with no hidden diagnostic fees.
- **Rechargeable Convenience:** (on applicable models) eliminates the need for constant battery replacements.

### Newman Medical: Your Trusted Partner in Ultrasound:

- **A Legacy of Innovation:** The Newman family has pioneered advancements in Doppler technology since the 1970s.
- **Founders of Industry Leaders:** Instrumental in the establishment of Imex Medical Systems (Nicolet/Natus) and Summit Doppler.
- **Commitment to Progress:** We continuously develop cutting-edge diagnostic solutions.
- **Your Success is Our Priority:** Focused on building long-term partnerships and exceeding expectations.

# DIGIDOP "CLASSIC" DD-700

HR-DISPLAY DOPPLER SYSTEM



## Accessories



Wall/Table Base  
(ACC-101)



Roll Stand  
(STND-101)



Carry Bag  
(ACC-170)



Ultrasound Gel  
8.8oz | 2.1oz  
(GEL-110) (GEL-100)

### ALWAYS BE READY

Ask about the rechargeable version of the DigiDop "Classic" HR-Display Fetal Doppler!



Smart Recharge System



Ask about our Doppler Recycling Program

## Ordering Information

### Doppler System

SKU	DESCRIPTION
DD-700	HR-Display Fetal Doppler System

### Probe Choices

SKU	DESCRIPTION
D2	2 MHz
D2W	2 MHz Waterproof
D3	3 MHz
D3W	3 MHz Waterproof
D5	5 MHz
D8	8 MHz



### Accessories

SKU	DESCRIPTION
ACC-101	Wall/Table Base
ACC-170	Carry bag for handheld Dopplers
STND-101	Roll stand with basket
GEL-110	Ultrasound gel 8.8oz 12/box
GEL-100	Ultrasound gel 2.1oz 12/box

## DD-700 Technical Specifications

<b>Dimensions:</b>	150 x 65 x 35 mm (6 x 2.5 x 1.25 in)
<b>Weight:</b>	340 grams (12 ounces)
<b>Battery type:</b>	3 x 1.5 nominal volts (AA)
<b>Level of protection against electrical shock:</b>	Type B Applied Part, Class II Equipment
<b>Comply with the following standards:</b>	IEC60101-1, IEC60601-2, IEC60601-2-37
<b>Operating temperature:</b>	10° ~ 40°C (50° ~ 104°F)
<b>Operating humidity:</b>	30% ~ 75%
<b>Transport/storage temperature:</b>	-20° ~ 50°C (-4° ~ 122°F)
<b>Transport/storage humidity:</b>	5% ~ 90%, non-condensing
<b>Battery life (AA):</b>	800 mins (Alkaline)
<b>Heart rate range:</b>	50 ~ 220 BPM
<b>Heart rate calculation accuracy:</b>	±3 BPM
<b>Sensitivity:</b>	9 weeks gestation (3MHz)
<b>Audio cable interface:</b>	3.5mm stereo plug
<b>Audio output power:</b>	1.1W
<b>GTIN:</b>	852785008883

## The DigiDop PRODUCT FAMILY

### DigiDop "Classic"



- Standard size
- Longer battery life than similar competitor products

### DigiDop II



- Smaller size
- Ergonomic - fits comfortably in your hand
- Similar battery life to competitor products

### DigiDop 901 Tabletop

- Ideal for hospital use
- Great room-filling sound from two speakers on top
- Two probe cradles
- Standard drop-on charging pyramid
- Optional wall mount recharger

