


**Fully assembled, mega racks are great for high throughput labs requiring large sample storage**

- Sturdy, durable material for easy and secure transport of clinical and research tubes
- Alphanumeric identification for quick location and retrieval of tubes
- Reusable, easy to clean
- Place different test groupings into one rack eliminating the need for multiple racks
- Stackable for staging and saving space
- Ideal footprint for use with programmable tube filling machines




120114

**Mega Racks® 10-13 mm - Holds 5-7 mL sample tubes**



| Item No. | Color   | Wells | Well Dia.<br>mm | L x W x H<br>in  | L x W x H<br>cm   | Material |
|----------|---------|-------|-----------------|------------------|-------------------|----------|
| HS24322A | Blue ■  | 216   | 13              | 13.3 x 9.3 x 3.6 | 33.7 x 23.5 x 9.1 | PP       |
| HS24322E | White □ |       |                 |                  |                   |          |
| 120112   | White □ | 432   | 13              | 26.5 x 9.3 x 3.5 | 67.3 x 23.5 x 9.0 | PP       |
| 120113   | Blue ■  |       |                 |                  |                   |          |



HS24323A & HS24323E

**Mega Racks® 13-16 mm - Holds 10 mL sample tubes**



| Item No. | Color   | Wells | Well Dia.<br>mm | L x W x H<br>in   | L x W x H<br>cm   | Material |
|----------|---------|-------|-----------------|-------------------|-------------------|----------|
| HS24323A | Blue ■  | 120   | 16              | 13.3 x 10.5 x 3.6 | 33.7 x 26.7 x 9.1 | PP       |
| HS24323E | White □ |       |                 |                   |                   |          |
| 120114   | White □ | 240   | 16              | 26.5 x 10.5 x 3.5 | 67.3 x 26.7 x 9.0 | PP       |
| 120115   | Blue ■  |       |                 |                   |                   |          |

**Mega Racks® 18-20 mm - Holds 15 mL conical/swab tubes and 20 mm round bottom tubes**



| Item No. | Color   | Wells | Well Dia.<br>mm | L x W x H<br>in   | L x W x H<br>cm   | Material |
|----------|---------|-------|-----------------|-------------------|-------------------|----------|
| 120834   | White □ | 120   | 20              | 14.9 x 11.9 x 3.6 | 37.7 x 30.1 x 9.1 | PP       |
| 120835   | Blue ■  |       |                 |                   |                   |          |



120834

**Mega Racks® 28-32 mm - Holds 50 mL conical/swab tubes**

| Item No. | Color   | Wells | Well Dia.<br>mm | L x W x H<br>in   | L x W x H<br>cm  | Material |
|----------|---------|-------|-----------------|-------------------|------------------|----------|
| *120958  | Blue ■  | 100   | 32              | 16.2 x 16.7 x 4.7 | 41.2 x 42.3 x 12 | ABS      |
| *120959  | White □ |       |                 |                   |                  |          |

\*Not recommended for autoclaving



120959