

ME2000

ELEVATING TREATMENT TABLE

STANDARD FEATURES:

- Load capacity: 400 lbs.
- Height range: 24-36"
- Steel base with white powder coat
- (2) Sliding doors on each side of table
- (1) Foot pedal
- Flat 2" high density foam cushion with nose hole
- Rounded corners



ME2000 with Lift Back Option



ME2000 with Lift Back and Split Leg Option

STANDARD TABLE:

ME2000

28"W x 77"L x 24-36"H

STANDARD UPHOLSTERY CHOICES:

	Royal Blue		American Beauty
	Imperial Blue		Lark
	Blue Ridge		Burgundy
	Deep Violet		Dove
	Forest		Charcoal
	Mandarin Orange		Mocha
	Black		Putty

OPTIONS:

- S2723** - 70° Lift Back Cushion
- S2590** - 45° Split Leg Cushion
- S2580** - Lift Back and Split Leg Cushion
- S2701** - Additional Foot Pedal
- S2779** - Additional Mounted Hand Control