

PTX HI-LO

ELECTRIC TREATMENT TABLE

STANDARD FEATURES:

- Load capacity: 550 lbs.
- 30"W x 77"L x 19-39"H
- Durable and stylish steel base (PTX9041)
- White anti-bacterial and anti-microbial powder coat on steel base
- Whisper quiet lift system with single foot control pedal
- Zero lateral displacement
- 30"W Cushions in 3 different configuration options
- 2" High density foam cushions with rounded corners
- Individual locking casters
- Nose hole
- Head section angle range: +40° to -90°
- Lower extremity angle range: 0 to 65°
- ADA Compliant



PTX1103

STANDARD TABLES:

PTX1101
1 Section
30" Width



PTX1102
2 Section
30" Width



PTX1103
3 Section
30" Width



OPTIONS:

- PT10001** - Nose Hole Plug
- S1077** - Dual Foot Control Pedal
- S4018** - Dual Foot Bar Control
- EPTPH01** - Paper Roll Holder

CUSTOM CUSHION BRANDING:

PT3030 - Logo Slip Cover
(30"W x 12"H)



STANDARD UPHOLSTERY CHOICES:

	Royal Blue		American Beauty
	Imperial Blue		Lark
	Blue Ridge		Burgundy
	Deep Violet		Dove
	Forest		Charcoal
	Mandarin Orange		Mocha
	Black		Putty