




American-made quality



Includes labor, parts, and accidental breakage



Includes cords and standard probes





Waterproof probes available



Rechargeable



Standard Options



Custom Option



Enhance your prenatal care with the DigiDop Fetal Doppler.

Designed for exceptional clarity and ease of use, the DigiDop empowers obstetricians to provide superior care throughout pregnancy. With advanced features and unmatched reliability, the DigiDop is the trusted choice for monitoring fetal heart rate and ensuring the well-being of both mother and baby.

Hear the Difference

- Digital Sound Clarity (DSC™) amplifies the signal up to 6x for exceptional accuracy.
- American-made with guaranteed quality.
- Durable and reliable with an unbreakable probe holder.
- Interchangeable probes for maximum efficiency.





Confidence-Inspiring Warranty

- 5-year coverage for all parts and labor, with no hidden fees.
- Rechargeable convenience (on applicable models).

Newman Medical: Your Ultrasound Partner

- Pioneers in Doppler technology since the 1970s.
- Founders of Imex Medical Systems and Summit Doppler.
- Committed to continuous development of cutting-edge solutions.
- Focused on building long-term partnerships and exceeding customer expectations.

DigiDop Professional and Sensitive Dopplers

-  Over **45 years in the marketplace**; a testament to our longevity and success
-  **Made in the USA** for superior quality and performance
-  **5-year all-inclusive warranty** for maximum protection and minimized downtime, headaches
-  **Best-in-class customer support**; we always answer the phone

DigiDop Classic

- Standard size
- Superior battery endurance compared to competitors



DigiDop II

- Compact size
- Comfortable palm-fitting design
- Equivalent battery endurance to competitors



DigiDop 901

- Ideal for hospital use
- Clear stereo sound from dual top speakers
- Two probe cradles
- Standard drop-on charging pyramid
- Optional wall mount recharger



Vascular Probes




Fetal Probes



Interchangeable probes can be used on any DigiDop

Vascular



- Great sensitivity
- Ergonomic shape easy to hold, place in sterile glove finger for use in sterile fields
 - **D8 – 8 MHz** - best for shallow vessels like radial, ankle, foot, digits, penile
 - **D5 – 5 MHz** - wider beam, deeper penetration. Good for brachial, carotid, femoral, deep veins

Fetal



- Great sensitivity
- Smooth shape, comfortable for patient, easy to hold and clean
 - **D3 – 3 MHz** - wide angle, ideal for early fetal hearts
 - **D3W** - a unique 3 MHz waterproof probe
 - **D2 – 2MHz** - deeper penetration
 - **D2W – 2 MHz** - waterproof

DPPG



- Unique audio PPG probe for use on digits. Great for toe-brachial studies

How to build your Doppler

DOPPLER CODE - PROBE CODE

	Doppler	Probe
DigiDop "Classic"	DD-300 Basic	D2 2MHz fetal probe
	DD-301 Basic + Recharge	D2W 2MHz waterproof fetal probe
	DD-700 HR Display	D3 3MHz fetal probe
	DD-701 HR Display + Recharge	D3W 3MHz waterproof fetal probe
	DD-901 Tabletop + HR Display + Recharge	D5 5MHz vascular probe
DigiDop II	DD-330 Basic	D8 8MHz vascular probe
	DD-330R Basic + Recharge	DPPG Audio PPG probe - includes toe cuff
	DD-770 HR Display	
	DD-770R HR Display + Recharge	