

PREMIER WOOD TREATMENT TABLE



STANDARD FEATURES:

- Load capacity: 550 lbs.
- 30"W x 78"L x 31"H
- Maple hardwood legs with plywood frame
- H-brace construction
- 2" High density foam cushion
- Rounded corners



APF-3078 with Nose Hole & Bottom Shelf

FLAT TABLES:

APF-3078 - 30" Flat Cushion
30"W x 78"L x 31"H

APFC-3078 - 30" Flat Cushion with Enclosed Cabinet
30"W x 78"L x 31"H

LIFT BACK TABLES:

APT-3078 - 30" Lift Back Cushion
30"W x 78"L x 31"H

APTC-3078 - 30" Lift Back Cushion with Enclosed Cabinet
30"W x 78"L x 31"H



APTC-3078

OPTIONS:

AP-Bottom Shelf - Bottom Shelf

AP-Cabinet - Shelf and Open Cabinet
Head End or Foot End

AP-NoseHole - Nose Hole in Cushion

AP-PaperCutter - Paper Roll Holder (On Foot End with Paper Cutter on Head End)

STANDARD UPHOLSTERY CHOICES:

	Royal Blue		American Beauty
	Imperial Blue		Lark
	Blue Ridge		Burgundy
	Deep Violet		Dove
	Forest		Charcoal
	Mandarin Orange		Mocha
	Black		Putty