

**Heart Rate Display, Alkaline Batteries, Pocket-Sized**  
with interchangeable probes for vascular and fetal applications

**INDUSTRY-LEADING  
5 YEAR WARRANTY**

**5** ALL PARTS & LABOR  
ACCIDENTAL BREAKAGE  
NO DIAGNOSTIC FEES  
NO HIDDEN FEES  
MINIMAL DOWN TIME

**MADE IN  
USA**

**REAL PEOPLE  
REAL SUPPORT**

NO ROBOTS.  
NO RUNAROUND.  
JUST REAL HELP.



The DigiDop II handheld Doppler combines advanced technology with ultimate portability. Made with pride in the USA using high-quality materials, this compact diagnostic tool delivers exceptional Digital Sound Clarity™ technology in a pocket-sized design. Its ergonomic form factor makes it perfect for mobile healthcare professionals.

Like all DigiDop products, it comes with an industry-leading 5-year warranty covering all parts and labor, with zero hidden fees or diagnostic charges.

## DigiDop II Features

- Made in the USA with high-quality materials
- Comprehensive 5-year warranty coverage
- Large, bright display
- Ultra-compact, pocket-sized design
- Ergonomic shape for comfortable handling
- Maximum portability without sacrificing performance
- Operates on 3 AAA alkaline batteries

## DD-770 Ordering Info

Model	Description	Pkg.
DD-770-D2	Heart Rate Display, Alkaline Batteries, 2 MHz fetal probe	1/ea
DD-770-D2W	Heart Rate Display, Alkaline Batteries, 2 MHz waterproof fetal probe	1/ea
DD-770-D3	Heart Rate Display, Alkaline Batteries, 3 MHz fetal probe	1/ea
DD-770-D3W	Heart Rate Display, Alkaline Batteries, 3 MHz waterproof fetal probe	1/ea
DD-770-D5	Heart Rate Display, Alkaline Batteries, 5 MHz vascular probe	1/ea
DD-770-D8	Heart Rate Display, Alkaline Batteries, 8 MHz vascular probe	1/ea

**Interchangeable probes**



Seeking a rechargeable alternative of the DD-770? The DD-770R delivers power on demand.

Smart Recharge System