


Disposable Cuvettes

- Fits most spectrophotometers and photometers
- Polystyrene for use in visible range (340-800nm)
- Protective packaging - dust-proof and non-scratch
- Absorption reproducibility error of less than $\pm 1\%$



MTC Bio's disposable cuvettes are compatible with most popular spectrophotometers. Manufactured from high quality, optically clear polystyrene, they are intended for determinations in the 340-800nm (visible) wavelength range. Available in three sizes: 1.5mL micro, 2.5mL semi-micro and 4.5mL macro. All sizes are 12.5mm square, 45mm height and are compatible with various spectrophotometer models including:

Carolina®, Spectronic®200, Spectronic®20 and 20D, Spectronic®Educator, Genesys™ and Accuris SmartDrop™.

Caps are available and can be ordered separately

SPECIFICATIONS:

Path Length:	10mm
Ext. Dimensions:	12.5 × 12.5 × 45mm
Wavelength Range:	340 to 800nm
Material:	Polystyrene
Temp range:	4°C to 70°C
Absorption:	Less than $\pm 1\%$

Round Opening Cuvettes



**1.5mL
S4100-1**



**2.5mL
S4100-2**

Square Opening Cuvettes



**1.5mL
S4100-3**



**4.5mL
S4100-4**

Disposable Cuvettes

Cat No.	Description	Qty
S4100-1	Cuvettes, PS, 1.5mL, Round Opening, 10 trays of 100 Cuvettes	1000 / pk
S4100-2	Cuvettes, PS, 2.5mL, Round Opening, 10 trays of 100 Cuvettes	1000 / pk
S4100-3	Cuvettes, PS, 1.5mL, Square Opening, 10 trays of 100 Cuvettes	1000 / pk
S4100-4	Cuvettes, PS, 4.5mL, Square Opening, 10 trays of 100 Cuvettes	1000 / pk

S4100-5: Round Opening Cap



S4100-6: Square Opening Cap



Cuvettes Caps

Cat No.	Description	Fits Cuvettes	Qty
S4100-5	Cuvette Plug Cap, Round Closure	<i>S4100-1, S4100-2</i>	500 / pk
S4100-6	Cuvette Plug Cap, Square Closure	<i>S4100-3, S4100-4</i>	500 / pk