



Cool Containers

Durable containers for chilling temperature-sensitive samples

Suitable for low temperature use. Ideal for ice, ice slurries and CO₂ pellets

- Made from polyurethane with metal inserts for stability
- Close-fit lid keeps evaporation to a minimum and features a comfortable recessed knob
- Tapered design for easy stacking
- Precise pour spout and large handles



Lid nests under bucket
for convenience and to save space

Item No.	Color	Style	Capacity liters	L x W x H in	L x W x H cm
HS28714K	Black	Bucket	4	7.1 (H) x 9.2 (D)	18.1 (H) x 23.4 (D)
HS28721G	Green	Pan	1	10.5 x 7.7 x 4.5	26.7 x 19.6 x 11.4
HS28724B	Blue	Pan	4	15 x 10 x 5.5	38.1 x 25.4 x 14
HS28729R	Red	Pan	9	20.5 x 11.5 x 6.4	52.1 x 29.2 x 16.3