

10 Year Limited
Warranty



Advantages

- ✓ **Powder Coated Fabrication**
Powder coating inside and out protects it from rust, extending the life of the scale
- ✓ **Exclusive 10 Year Limited Warranty**
Health o meter® Professional's industry leading warranty protects you against defects of materials or workmanship



Model	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
400KL	58.875" x 20.25" x 6.375" / 4.40	28.75	892076002121

Specifications

- Capacity: 390 lb / 180 kg
- Graduation: 1/4 lb / 100 g
- Platform Size: 10-1/2" (w) x 14" (d) x 3-1/4" (h), 267 mm (w) x 356 mm (d) x 83 mm (h)
- Optional height rod available (item # 402KLROD)
- Optional wheels available (item # 55000)
- Optional counter weights available, to increase capacity to 490lb / 210kg (item # 55090)
- 10 Year Limited Warranty

Optional Wheels
Item # 55000

