

# Health o meter®

P r o f e s s i o n a l

10 Year Limited  
Warranty



## Advantages

- ✓ **Powder Coated Fabrication**  
Powder coating inside and out protects it from rust, extending the life of the scale
- ✓ **Height Rod Included**  
The easy-to-use height rod provides height measurements from 11 7/8" - 84"/30cm - 213cm
- ✓ **Exclusive 10 Year Limited Warranty**  
Health o meter® Professional's industry leading warranty protects you against defects of materials or workmanship

**Health o meter®**  
P r o f e s s i o n a l  
[www.homscales.com](http://www.homscales.com)

**Health o meter®**  
P r o f e s s i o n a l



| Model | Shipping Dims. (L x W x H) / Cube | Weight (LB) | UPC Code     |
|-------|-----------------------------------|-------------|--------------|
| 402KL | 58.875" x 20.25" x 6.375" / 4.40  | 33.25 lb    | 892076002015 |

**Optional Wheels**  
Item # 55000



**Optional Counterweights**  
Item # 55090



**Specifications**

- Capacity: 390 lb / 180 kg
- Graduation: 1/4 lb / 100 g
- Platform Size (w x d x h): 10 1/2" x 14" x 3 1/4" / 267 mm x 356 mm x 83 mm, platform mat: non-slip black
- Fixed Beam, pounds and kilograms can be read on both sides
- Professional Dual Reading Height Rod: 11 7/8" - 84"/30cm - 213cm
- Optional wheels available (item # 55000)
- Optional counter weights available, to increase capacity to 490lb / 210kg (item # 55090)
- 10 Year Limited Warranty