


Health o meter®
P r o f e s s i o n a l

10 Year Limited
Warranty



Advantages

- ✓ **Adjustable Unit Measurement Bar**
The mechanical rotating poise bars allow for easy reading of weight in lb or kg
- ✓ **Exclusive 10 Year Limited Warranty**
Health o meter® Professional's industry leading warranty protects you against defects of materials or workmanship
- ✓ **Height Rod Included**
The easy-to-use height rod provides height measurements from 23" to 84"
- ✓ **Powder Coated Fabrication**
Powder coating inside and out protects it from rust, extending the life of the scale
- ✓ **Engraved Numerals**
Engraved numerals on the poise bars can withstand frequent use, providing long-lasting readability

Health o meter® P r o f e s s i o n a l



Exclusive Engraved Numerals

Easy to read rotating bars



Powder coated inside and out for extended life

Easy to read rotating bars



Exclusive rotating poise bar

Model	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
450KL	58.875" x 20.25" x 6.375" / 4.40	36.75 lb	892076002169

Specifications

- Capacity: 500 lb / 200 kg – with included counterweights.
- Graduation: 1/4 lb / 100 g
- Platform Size: 10 1/2" (w) x 14" (d) x 3 1/4" (h) / 267 mm (w) x 356 mm (d) x 83 mm (h)
- Height Rod: 23 5/8"- 84" / 60 cm - 213 cm
- Wheels: Optional item #55000
- 10 Year Limited Warranty

Optional Wheels Item # 55000

