



This classic “Bigfoot” scale has a deep platform with non-slip surface and an 8” easy-to-read raised dial—good for patients with vision challenges or those who prefer more foot room.

Advantages

- ✓ **180 kg Capacity Floor scale**
The high capacity of this floor scale provides more versatile use
- ✓ **8" Raised Dial Display**
Large, easy-to-read dial display is raised above the scale base to allow for more footroom on the platform
- ✓ **Weighs in KG**
Dial display shows weight in kilograms



Specifications

- Capacity: 180 kg
- Graduation: 1 kg
- 8" dial (raised 4")
- Platform Size (w x d x h):
11" x 12 ½" x 3" / 279 mm x 318 mm x 76 mm
- Product Footprint (w x d x h):
11 ¾" x 18 ¼" x 7 ¼" / 298 mm x 464 mm x 184 mm
- Product Weight: 9 lb / 4 kg
- 1 Year Limited Warranty

Model	Description
160KL	Mechanical floor scale - kilograms & pounds (2 pack)
160LB	Mechanical floor scale - pounds only (2 pack)
160KG	Mechanical floor scale - kilograms only (2 pack)
160KLS	Mechanical floor scale - kilograms & pounds (single)
160LBS	Mechanical floor scale - pounds only (single)
160KGS	Mechanical floor scale - kilograms only (single)



This classic “Bigfoot” scale has a deep platform with non-slip surface and an 8” easy-to-read raised dial—good for patients with vision challenges or those who prefer more foot room.

Advantages

- ✓ **400 lb / 180 kg Capacity Floor scale**
The high capacity of this floor scale provides more versatile use
- ✓ **8" Raised Dial Display**
Large, easy-to-read dial display is raised above the scale base to allow for more footroom on the platform
- ✓ **Weighs in LB/KG**
Dial display shows weight in both pounds and kilograms providing versatility

Health o meter® P r o f e s s i o n a l



Specifications

- Capacity: 400 lb / 180 kg
- Graduation: 1 lb / 1 kg
- 8" dial (raised 4")
- Platform Size (w x d x h):
11" x 12 ½" x 3" / 279 mm x 318 mm x 76 mm
- Product Footprint (w x d x h):
11 ¾" x 18 ¼" x 7 ¼" / 298 mm x 464 mm x 184 mm
- Product Weight: 9 lb / 4 kg
- 1 year limited warranty



This classic “Bigfoot” scale has a deep platform with non-slip surface and an 8” easy-to-read raised dial—good for patients with vision challenges or those who prefer more foot room.

Advantages

- ✓ **400 lb Capacity Floor scale**
The high capacity of this floor scale provides more versatile use
- ✓ **8" Raised Dial Display**
Large, easy-to-read dial display is raised above the scale base to allow for more footroom on the platform
- ✓ **Weighs in LB**
Dial display shows weight in pounds

Health o meter® P r o f e s s i o n a l



Specifications

- Capacity: 400 lb
- Graduation: 1 lb
- 8" dial (raised 4")
- Platform Size (w x d x h):
11" x 12 ½" x 3" / 279 mm x 318 mm x 76 mm
- Product Footprint (w x d x h):
11 ¾" x 18 ¼" x 7 ¼" / 298 mm x 464 mm x 184 mm
- Product Weight: 9 lb / 4 kg
- 1 Year Limited Warranty