

Capacity: 330 lb



Advantages

- 330 lb Capacity Floor scale**
 The high capacity of this floor scale provides more versatile use
- Easy-to-read Dial**
 Large, easy-to-read dial is conveniently located at the top of the scale
- Durable Construction**
 Solid construction and durable, non-skid mat can withstand frequent use
- Small-Step Platform**
 Low platform height requires only minimal step up, offering patient a more comfortable and stable weighing experience

Capacity: 330 lb



175LB Dial

Specifications

- Capacity: 330 lb
- Graduation: 1 lb
- Platform Size (w x d x h):
10 ½" x 9 ⅜" x 2" / 266 mm x 238 mm x 51 mm
- Product Footprint (w x d x h):
11 ⅜" x 9 ¾" x 2" / 289 mm x 248 mm x 57 mm
- Optional Carrying Case Available
(order # 64771)
- 1 Year Limited Warranty



Carrying Case for 175LB
Item# 64771