

Health o meter[®]

P r o f e s s i o n a l

EXCLUSIVE ANTI-SLIP TREAD!



Advantages

- Exclusive Anti-Slip Tread**
 Chain-link bonded surface tread provides additional traction and peace of mind
- Medical Grade Home Care Scale**
 Designed for use in the home/bathroom, with accuracy and features of a professional scale
- 440 lb / 200 kg Capacity Floor Scale**
 The high capacity of this floor scale provides more versatile use
- Backlit Display**
 Easy to read from any angle and in any light
- 5-Second Weight Hold**
 Weight measurement will stay displayed for 5 seconds after stepping off the scale allowing the user or caregiver time to stabilize themselves or the patient

Health o meter[®]
 P r o f e s s i o n a l
www.homscales.com

See reverse side for additional specifications

Health o meter®

P r o f e s s i o n a l



Specifications

- Capacity: 440 lb / 200 kg
- Resolution: 0.2 lb / 0.1 kg
- Platform Size (w x d x h):
13" x 13" x 7/8" / 330 mm x 330 mm x 22 mm
- Product Footprint (w x d x h):
13" x 13" x 7/8" / 330 mm x 330 mm x 22 mm
- Product Weight: 4 lb / 1.8 kg
- Display: 1 7/8" / LCD
- Power Source: (2) 3-volt lithium batteries (included)
- Functions: LB/KG Button, 5 Second Weight Hold, Auto Off, Auto Zero
- 1 year limited warranty