



### Advantages

- ✓ **400 lb/180 kg Capacity Floor Scale**  
The high capacity of this floor scale provides more versatile use
- ✓ **1 ½" LCD Display**  
The large display shows the user an easy-to-read weight measurement
- ✓ **5-Second Weight Hold**  
Weight measurement will stay displayed for 5 seconds after stepping off scale allowing the user time to stabilize themselves or the patient



Carrying Case for 800KL  
Item# 64771



### Specifications

- Capacity: 400 lb / 180 kg
- Resolution: 0.1 lb / 0.05 kg
- 1 ½" LCD display
- Functions: LB / KG Switch, Auto Zero, 5 Second Weight Hold, Auto Off
- Platform Size (w x d x h):  
13 ¾" x 10 ¾" x 1 ½"  
349 mm x 273 mm x 38 mm
- Product Footprint (w x d x h):  
14 ¼" x 11 ⅞" x 1 ½"  
362 mm x 302 mm x 38 mm
- (2) AAA batteries included
- Optional Carrying Case Available (order # 64771)
- 1 Year Limited Warranty