



This scale can weigh immobile patients up to 400 lb / 181 kg. It has a durable steel structure, fabric sling, and wheels for portability.

Advantages

- EMR Connectivity**
Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter® Professional Connectivity Solutions
- Value-Added Features**
Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Pre-tare Hold/Release, Reweigh, Recall, Auto Zero, Auto Off) included at no additional cost
- Antimicrobial Keypad**
The keypad features antimicrobial protection
- Calculates Body Mass Index (BMI)**
Scale functions include the ability to calculate the patient's BMI
- Hospital Grade Construction**
Durable steel structure with large sturdy stretcher enables patients to be weighed accurately and comfortably
- Everlock® (patent pending)**
Option that allows users to permanently lock the scale in either LB or KG to improve safety

Healthometer®

P r o f e s s i o n a l



Specifications



- Capacity: 400 lb / 181 kg
- Resolution: 0.2 lb / 0.1 kg
- Connectivity via USB, optional Pelstar® wireless technology
- Display: 1 ½" / 38 mm High-Contrast Color TFT-LCD Screen
- Canvas Size (w x d):
70 ½" x 30"
1791 mm X 762 mm
- Product Footprint - at highest extension of lift (w x d x h):
70" x 44 ½" x 77"
1778 mm x 1130 mm x 1956 mm
- Product Weight: 104 lb / 47 kg
- Functions: LB / KG Conversion, LB / KG lock Out, BMI, Zero, Tare, Pre-tare, Hold / Release, Reweigh, Auto Zero, Auto Off, Everlock®
- 4 Wheels for ease of movement, 1 locking wheel for safety
- Power Source: 100-240V adapter (ADPT30) included, 6 D batteries (not included)
- 2 year limited warranty

Model	Description
2000KL	Digital in-bed / stretcher scale, pounds and kilograms
2000KL-BT	2000KL with Pelstar wireless technology
2000KG	Digital in-bed / stretcher scale, kilograms only
2000KG-BT	2000KG with Pelstar wireless technology
SS-2000KL	ScaleSurance 2 year extended warranty for 2000KL, 2000KG, 2000KL-BT & 2000KG-BT*
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
PELSTARONE-P	Wireless thermal printer