



3400KL (pan)



3401KL (tray)

The solid, industrial design of the 3400KL includes a stable, stainless steel base and removable, stainless steel weighing pan (3401KL comes with a flat tray) that can be autoclaved for sterilization.

## Advantages

- Stainless Steel Base and Weighing Surface**  
The stainless steel design is easy to clean and ideal for sterile environments.
- Rechargeable Battery Power Source**  
Scale is powered by a 6V rechargeable battery with battery charger included with the scale.
- High Resolution Measurements**  
Scale provides highly precise and accurate results with refined 0.05 oz / 1 g resolution.
- Backlit Display**  
Easy to read from any angle and in any light condition.
- Convenient Weighing Pan (3400KL)**  
The 2½" deep and 9½" wide stainless steel weighing pan effectively retains fluids and can be removed for sterilization.
- Convenient Weighing Tray (3401KL)**  
9¾" x 7¾" tray is useful for weighing items that require a flat surface. Stainless steel tray is removable and can be autoclaved.
- Tare Functionality**  
Container weight can be easily deducted to provide accurate net weight of the contents.



Optional additional pans and trays are available, for pans order 3400PAN and for trays order 3400TRAY

**Specifications**

- Capacity: 176 oz / 5000 g
- Resolution: 0.05 oz / 1g
- Weighing Surface Dimensions (pan) (w x d x h): 9 1/2" x 11 3/4" x 2 1/2" / 241mm x 298mm x 64mm
- Weighing Surface Dimensions (tray) (w x d x h): 7 3/4" x 9 3/4" x 5/8" / 197mm x 248mm x 16mm
- Scale Base Dimensions (including handles) (w x d x h): 11 1/2" x 9 7/8" x 5 1/4" / 292mm x 251mm x 133mm
- Product Footprint 3400KL (w x d x h): 10 7/8" x 13" x 8" / 276mm x 330mm x 203mm
- Product Footprint 3401KL (w x d x h): 10" x 11 3/4" x 5 3/8" / 254mm x 298mm x 137mm
- Product Weight 3400KL: 11 lb / 5 kg
- Product Weight 3401KL: 10 lb / 4.5 kg
- Display: 5/8" (16mm) LCD Display
- Power Source: 100 – 240 VAC battery charger / 6 volt rechargeable battery (included)
- Functions: Zero Out/Tare, Auto Zero, Auto Off, g/oz Conversion, Low Battery Indicator, Backlit Display
- 2 Year Limited Warranty

**Options & Accessories**

Model	Description
3400KL	Digital wet diaper / lap sponge / organ scale with pan
3401KL	Digital wet diaper / lap sponge / organ scale with tray
SS-3400KL	ScaleSurance 2 year extended warranty for 3400KL
SS-3401KL	ScaleSurance 2 year extended warranty for 3401KL
3400PAN	3 pack of stainless steel pans for 3400KL & 3401KL
3400TRAY	3 pack of stainless steel trays for 3400KL & 3401KL
3400ADPT	Battery charger for 3400KL & 3401KL