



Issuing Date: 7-Jan-1998

Revision Date: 5-Jun-2018

Rev. 5

1. IDENTIFICATION OF THE SUBSTANCE/ MIXTURE AND THE COMPANY/UNDERTAKING

Product Name: Kanamycin Sulfate
Product Code: IB02120
Recommend Use: For further manufacturing use only
Uses Advised Against: Not for human or animal drug use
Synonyms: CANTREX * CRISTALOMICINA * KAMYCIN * KAMYNEX * KAMANABRISTOL * KANACEDIN * KANAMYCIN A SULFATE * KANAMYCIN MONOSULFATE * KANAMYCINSULFATE * KANAMYTREX * KANAQUA * KANASIG * KANATROL * KANESCIN * KANICIN * KANNASYN * KANO * KANTREX * KANTREXIL * KANTROX * KLEBCIL * OPHTALMOKALIXAN * OTOKALIXIN * RESISTOMYCIN (BAYER) *
Company: IBI SCIENTIFIC
7445 Chavenelle Road
Dubuque, IA 52002
Emergency Tele. Number: 3E Company (Acct# 9027) 1-800-451-8346 / (001) 760-602-8703
Company Tele. Number: 1-800-253-4942 / (001) 563-690-0484
E-Mail Address: info@ibisci.com

2. HAZARDOUS IDENTIFICATION

GHS – Classification

OSHA Regulatory Status: This chemical is considered hazardous by the 2012 OSHA Hazard Communication Standard (29CFR 1910.1200)

Reproductive Toxicity: Category 1B

Signal Word: Danger



Appearance: Fine white to off-white

Physical State: Powder

Odor: No information available

Hazard Statements: H360-May damage fertility or the unborn child

Precautionary Statements: P201-Obtain special instructions before use
P202-Do not handle until all safety precautions have been read and understood
P281-Use personal protective equipment as required
P308+P313-If exposed or concerned: Get medical advice/attention
P405-Store locked up
P501-Dispose of contents/container to an approved waste disposal plant

Hazards not otherwise classified: Not applicable

Other Information: May be harmful if swallowed

3. COMPOSITION/INFORMATION ON INGREDIENTS

Chemical Name: Kanamycin Sulfate
EC No.: EEC No. 246-933-9
REACH Reg. No.: Not Available
CAS No.: 25389-94-0
Weight %: 95 – 100
Classification: R43

4. FIRST AID MEASURES

Eye Contact: Rinse thoroughly with plenty of water for at least 15 minutes, lifting lower and upper eyelids. Consult a physician.
Skin Contact: Remove contaminated clothing and shoes. In case of contact with substance, wipe from skin immediately, flush skin or eyes with running water for at least 20 minutes.
Inhalation: Move to fresh air. If breathing becomes difficult, call a physician.
Ingestion: Clean mouth with water and afterwards drink plenty of water. Never give anything by mouth to an unconscious person.
Notes to Physician: Treat symptomatically.

5. FIRE FIGHTING MEASURES

Flammable Properties: Not flammable
Flash Point: Not determined
Suitable Extinguishing Media: Dry chemical, CO₂, water spray, regular foam
Hazardous Combustion Products: Carbon oxides, Nitrogen oxides (NO_x),
Explosion Data
Sensitivity to Mechanical Impact: Not sensitive
Sensitivity to Static Discharge: Not sensitive
Protective Equip. & Precautions for Firefighters: As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equiv.) and full protective gear.

6. ACCIDENTAL RELEASE MEASURES

Personal Precautions: Avoid contact with skin, eyes, and clothing. Ensure adequate ventilation.
Environmental Precautions: See Section 12 for additional ecological information.
Methods for Containment: Prevent further leakage or spillage, if safe to do so.
Methods for Cleaning Up: Avoid dust formation. Pick up and transfer to properly labeled containers. Ventilate area and wash spill site after material pickup is complete.

7. HANDLING AND STORAGE

Handling: Handle in accordance with good industrial hygiene and safety practice.
Storage: Keep container tightly closed in a dry, cool, and well-ventilated place.
Incompatible Products: None known based on supplied information.

8. EXPOSURE CONTROLS / PERSONAL PROTECTION

Exposure Guidelines: This product does not contain any hazardous materials with occupational exposure limits established by the region specific regulatory bodies.

Chemical Name: Kanamycin sulfate - #25389-94-0

ACGIH TLV: -

OSHA PEL: -

NIOSH IDLH: -

Engineering Measures: Showers, Eyewash Stations, Ventilation Systems

Personal Protective Equipment

Eye/Face Protection: Safety glasses with side-shields.

Skin & Body Protection: Wear protective gloves/clothing.

Respiratory Protection: If exposure limits are exceeded or irritation is experienced, NIOSH/MSHA approved respiratory protection should be worn. Positive-pressure supplied air respirators may be required for high airborne contaminant concentrations. Respiratory protection must be provided in accordance with current local regulations.

Hygiene Measures: Handle in accordance with good industrial hygiene and safety practice.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance: Fine white to off-white
Odor: No information available
Odor Threshold: No information available
Decomposition Temp. °C: No information available
Freezing Point: No information available
Initial Boiling Point: No information available
Physical State: Powder
pH: No data available
Flash Point: Not determined
Autoignition Temp. °C: No data available
Boiling Point/Range: No data available
Melting Point/Range: No information available
Flammability Limits (Air): Upper-No data available * Lower-No data available
Explosion Limits No information available
Explosive Properties: No information available
Oxidizing Properties: No information available
Evaporation Rate: No information available
MMHG @ 37.8°C: No data available
Vapor Density: No data available
Specific Gravity: No data available
Water Solubility: No data available
Solubility in other solvents: No data available
Partition Coefficient (n-octanol/water): No data available
Viscosity: No information available

10. STABILITY AND REACTIVITY

Reactivity:	No data available
Chemical Stability:	Stable under recommended storage conditions
Possibility of Hazardous Reactions:	None under normal processing
Incompatible Products:	None known based on information supplied
Conditions to Avoid:	Extremes of temperature and direct sunlight.
Hazardous Decomposition Products:	Carbon oxides, Nitrogen oxides (NOx)
Hazardous Polymerization:	Hazardous polymerization does not occur.

11. TOXICOLOGICAL INFORMATION

Information on likely routes of exposure

Product Information:

Inhalation:	There is no data available for this product
Eye Contact:	There is no data available for this product
Skin Contact:	There is no data available for this product
Ingestion:	There is no data available for this product

Chemical Name:	Kanamycin sulfate - #25389-94-0
LD50 Oral:	>4mg/kg (Rat)
LD50 Dermal:	-
LC50 Inhalation:	-

Information on toxicological effects

Symptoms: No information available.

Delayed and immediate effects as well as chronic effects from short and long term exposure

Sensitization:	No information available
Mutagenic Effects:	No information available
Carcinogenicity:	No information available
Chemical Name:	Kanamycin sulfate - #25389-94-0
ACGIH:	-
IARC:	-
NTP:	-
OSHA:	-

Reproductive Toxicity: May cause harm to the unborn child

STOT – single exposure: No information available

STOT – repeated exposure: No information available

Aspiration Effects: No information available

Numerical measures of toxicity – Product Information

The following values are calculated based on chapter 3.1 of the GHS document

12. ECOLOGICAL INFORMATION

Ecotoxicity:	100% of mixture consists of component(s) of unknown hazards to aquatic environment
Chemical Name:	Kanamycin sulfate - #25389-94-0
Toxicity to Algae:	-
Toxicity to Fish:	-
Daphnia Magna (Water Flea):	-
Persistence and Degradability:	No information available
Bioaccumulation/Accumulation:	No information available
Chemical Name:	Kanamycin sulfate - #25389-94-0
Low Pow:	-
Other Adverse Effects:	No information available

13. DISPOSAL CONSIDERATIONS

Water Disposal Method: Dispose of material in accordance with all federal, state, and local regulations.

Contaminated Packaging: Do not re-use empty containers.

Chemical Name: Kanamycin sulfate - #25389-94-0

RCRA – Halogenated Organic Compounds: -

RCRA – P Series Wastes: -

RCRA – F Series Wastes: -

RCRA – K Series Wastes: -

Chemical Name: Kanamycin sulfate - #25389-94-0

California Hazardous Waste Status: -

14. TRANSPORT INFORMATION

IATA
Not regulated

US DOT
Not regulated

15. TOXICOLOGICAL INFORMATION

International Inventories

TSCA: Does not comply

DSL/NDSL: Does not comply

EINECS/ELINCS: Complies

ENCS: Does not comply

IECSC: Complies

KECL: Complies

PICCS: Complies

AICS: Complies

Legend:

TSCA: United States Toxic Substances Control Act Section 8(b) Inventory

DSL/NDSL: Canadian Domestic Substances List/Non-Domestic Substances List

EINECS/ELINCS: European Inventory of Existing Commercial Chemical Substances/EU List of Notified Chemical Substances

ENCS: Japan Existing and New Chemical Substances

IECSC: China Inventory of Existing Chemical Substances

KECL: Korean Existing and Evaluated Chemical Substances

PICCS: Philippines Inventory of Chemicals and Chemical Substances

AICS: Australian Inventory of Chemical Substances

U.S. Federal Regulations

SARA 313 Section 313 of Title III of the superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain a chemical or chemicals which are subject to the reporting requirements of the Act & Title 40 of the Code of Federal Regulations, part 372.

SARA 311/312 Hazard Categories

Acute Health Hazard No

Chronic Health Hazard Yes

Fire Hazard No

Sudden Release of Pressure Hazard No

Reactive Hazard No

Clean Water Act This product does not contain any substances regulated as pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42).

Chemical Name: Kanamycin sulfate - #25389-94-0
CWA – Reportable Quantities: -
CWA – Toxic Pollutants: -
CWA – Priority Pollutants: -
CWA – Hazardous Substances: -

CERCLA This material, as supplied, does not contain any substances regulated as hazardous substances under the Comprehensive Environmental Response Compensation and Liability Act (CERCLA) (40CFR 302) or the Superfund Amendments and Reauthorization Act (SARA) (40CFR 355). There may be specific reporting requirements at the local, regional, or state level pertaining to releases of this material.

Chemical Name: Kanamycin sulfate - #25389-94-0
Hazardous Substances RQs: -
Extremely Hazardous Substances RQs: -
RQ: -

U.S. State Regulations

California Proposition 65 This product does not contain any Proposition 65 chemicals.

Chemical Name: Kanamycin sulfate - #25389-94-0
California Prop. 65: -

U.S. Right-to-Know Regulations

Chemical Name: Kanamycin sulfate - #25389-94-0
New Jersey: -
Massachusetts: -
Pennsylvania: -

U.S. EPA Label Information

EPA Pesticide Registration #: Not applicable.

Canada No information available

WHMIS Hazard Class

D2A – Very toxic materials



16. OTHER INFORMATION

Issuing Date:	7-Jan-1997
Revision Date:	5-Jun-2018
Revision Note:	Updated address and logo
Recommended Restrictions:	No information available
Disclaimer:	The information provided on this SDS is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guide for safe handling, use, processing, storage, transportation, disposal, and release; and is not to be considered as a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other material or in any process, unless specified in the test.

End of Safety Data Sheet