

ESSENTIAL THERA-P™

ELECTRIC TREATMENT TABLE

STANDARD FEATURES:

- Load capacity: 650 lbs.
- 27"W x 77"L x 17-35"H
- Durable and stylish steel base (EPT9020/EPT9028)
- White anti-bacterial and anti-microbial powder coat base
- Whisper quiet lift system with single foot control pedal
- 27"W Cushions in 7 different configuration options
- 2" High density foam cushions with rounded corners
- Individual locking casters
- Head section angle range: +40° to -90° (EPT77: +35° to -45°)
- Lower extremity angle range: 0 to 65° (EPT77: 0 to 75°)
- Armrest horizontal range: 0 to 45°
- Armrest vertical range: 0 to 7"
- 7 Section configuration includes electric elevating midsection with hand control
- ADA Compliant




EPT03

NOTE: International (240V) option available

TABLE CONFIGURATIONS:



EPT01
1 Section
27" Width



EPT02
2 Section
27" Width



EPT03
3 Section
27" Width




EPT04
4 Section
27" Width



EPT05
5 Section
27" Width



EPT06
6 Section
27" Width



EPT77
7 Section
27" Width

TABLE CONFIGURATIONS WITH ELECTRIC LIFT BACK

EPT02-EL
2 Section
27" Width

EPT03-EL
3 Section
27" Width

EPT05-EL
5 Section
27" Width

TABLE FEATURES:

<p>Flat Cushion</p> <p>EPT01</p>	<p>Lower Extremity Angle Ranges from 0 to 65°</p> <p>EPT02/EPT02-EL</p>	<p>Lower Extremity Angle Ranges from 0 to 65°</p> <p>Optional Elevating Midsection</p> <p>Head Section Angle Ranges from +40° to -90°</p> <p>EPT03/EPT03-EL</p>	<p>Lower Extremity Angle Ranges from 0 to 65°</p> <p>Head Section Angle Ranges from +40° to -90°</p> <p>EPT04</p>
<p>Lower Extremity Angle Ranges from 0 to 65°</p> <p>Each Armrest Vertical Ranges from 0 to 7"</p> <p>Optional Elevating Midsection</p> <p>Head Section Angle Ranges from +40° to -90°</p> <p>EPT05/EPT05-EL</p>	<p>Lower Extremity Angle Ranges from 0 to 65°</p> <p>Each Armrest Vertical Ranges from 0 to 7"</p> <p>Head Section Angle Ranges from +40° to -90°</p> <p>Each Armrest Horizontal Ranges from 0 to 45°</p> <p>EPT06</p>	<p>Lower Extremity Angle Ranges from 0 to 75°</p> <p>Standard Electric Elevating Midsection</p> <p>Shoulder Dropouts 0 to -90°</p> <p>Head Section Angle Ranges from +35° to -45°</p> <p>Adjustable Armrests</p> <p>EPT77</p>	

OPTIONS:

OPTIONS:

EPT9021 - 650 lb Lifting Capacity (240V FOR EPT01, EPT02, EPT03, EPT04, EPT05, EPT06)

EPT9029 - 650 lb Lifting Capacity (240V FOR EPT77)

EPT303X - 30"W Cushion Configuration (not available on EPT77)

S1135 - Manual Elevating Midsection (only available on EPT03 & EPT05)

S1077 - Dual Foot Control Pedal

S1078 - Dual Foot Bar Control

S2078 - Dual Foot Bar Control With Front Foot Bar Adder

S1151 - (1) Battery & (1) Wall Charger

S1152 - (2) Batteries & (1) Wall Charger

EPTPH01 - Paper Roll Holder (only available on EPT01, EPT02, or EPT03)


PT10001 - Nose Hole Plug



Manual Elevating Midsection (S1135)



Dual Foot Control Pedal (x2) (S1077)



Dual Foot Bar Control (S1078)



Dual Foot Bar Pedal With Front Foot Bar Adder (S2078)




(1) Battery & (1) Wall Charger (S1151)



(2) Batteries & (1) Wall Charger (S1152)



Paper Roll Holder (EPTPH01)



Nose Hole Plug (PT10001)

STANDARD FRAME COLORS:




White

CUSTOM BASE COLORS:



STANDARD UPHOLSTERY CHOICES:



CUSTOM BASE OPTIONS:

S1178 - Table Base Graphics

PT9020-CC- Pick Custom Base Color

S2000 - Table Base Powder Coat Color and Table Base Graphics

