

THERA-P™

ELECTRIC TREATMENT TABLE

STANDARD FEATURES:

- Load capacity: 650 lbs.
- 27"W x 77"L x 17-35"H
- Durable and stylish steel base (PT9020/PT9028)
- White anti-bacterial and anti-microbial powder coat base
- Whisper quiet lift system with single foot control pedal
- Height adjustment: 17-35"
- 27"W Cushions in 8 different configuration options
- 2" High density foam cushions with rounded corners
- Low profile locking casters
- Head section angle range: +40° to -90° (PT1177: +35° to -45°)
- Lower extremity angle range: 0 to 65° (PT1177: 0 to 75°)
- Armrest horizontal range: 0 to 45°
- Armrest vertical range: 0 to 7"
- 7 Section configuration includes electric elevating midsection with hand control
- ADA Compliant



PT1104-SL with Dual Foot Bar
Control and Branding Kit 3

NOTE: International (240V) option available

TABLE CONFIGURATIONS WITH MANUAL LIFT BACK:




PT1101
1 Section
27" Width



PT1102
2 Section
27" Width




PT1103
3 Section
27" Width




PT1104
4 Section
27" Width



PT1105
5 Section
27" Width




PT1106
6 Section
27" Width



PT1177
7 Section
27" Width

PT1104-SL
4
Section with
Split Leg 30"
Width

TABLE CONFIGURATIONS WITH ELECTRIC LIFT BACK:



PT1102-EL
2 Section
27" Width



PT1103-EL
3 Section
27" Width



PT1105-EL
5 Section
27" Width

STANDARD UPHOLSTERY CHOICES:

	Royal Blue		American Beauty
	Imperial Blue		Lark
	Blue Ridge		Burgundy
	Deep Violet		Dove
	Forest		Charcoal
	Mandarin Orange		Mocha
	Black		Putty







OPTIONS:

TABLE FEATURES:

 <p>Flat Cushion</p>	 <p>Lower Extremity Angle Ranges from 0 to 65°</p>	 <p>Lower Extremity Angle Ranges from 0 to 65° Optional Elevating Midsection Head Section Angle Ranges from +40° to -90°</p>	 <p>Lower Extremity Angle Ranges from 0 to 65° Head Section Angle Ranges from +40° to -90°</p>
PT1101	PT1102/PT1102-EL	PT1103/PT1103-EL	PT1104
 <p>*Table shown with options</p>	 <p>Lower Extremity Angle Ranges from 0 to 65° Each Armrest Vertical Ranges from 0 to 7" Head Section Angle Ranges from +40° to -90° Optional Elevating Midsection Each Armrest Horizontal Ranges from 0 to 45°</p>	 <p>Lower Extremity Angle Ranges from 0 to 65° Each Armrest Vertical Ranges from 0 to 7" Head Section Angle Ranges from +40° to -90° Each Armrest Horizontal Ranges from 0 to 45°</p>	 <p>Lower Extremity Angle Ranges from 0 to 75° Standard Electric Elevating Midsection Head Section Angle Ranges from +35° to -45° Adjustable Armrests Shoulder Dropouts 0 to -90°</p>
PT1104-SL	PT1105/PT1105-EL	PT1106	PT1177

OPTIONS:

PT9021 - 650 lb Lifting Capacity (240V FOR PT1101, PT1102, PT1103, PT1104, PT1105, PT1106)
PT9029 - 650 lb Lifting Capacity (240V FOR PT1177)
PT303X - 30"W Cushion Configuration (not available on PT1177)
S1135 - Manual Elevating Midsection (only available on PT1103 & PT1105)
S1077 - Dual Foot Control Pedal
S1078 - Dual Foot Bar Control
S2078 - Dual Foot Bar Control With Front Foot Bar Adder
S1151 - (1) Battery & (1) Wall Charger
S1152 - (2) Batteries & (1) Wall Charger
EPTPH01 - Paper Roll Holder (only available on PT1101, PT1102, or PT1103)
PT1001 - Nose Hole Plug

		
Manual Elevating Midsection (S1135)	Dual Foot Control Pedal (x2) (S1077)	Dual Foot Bar Control (S1078)
		
Dual Foot Bar Pedal With Front Foot Bar Adder (S2078)	(1) Battery & (1) Wall Charger (S1151)	Paper Roll Holder (EPTPH01)

BRANDING KITS:

BRANDING KIT 1 - S1170 · Slip Cover · Custom Table Base Color · Custom Table Base Graphics	BRANDING KIT 4 - S1174 · Deboss Logo · White Table Base Color · Custom Table Base Graphics
BRANDING KIT 2 - S1171 · Slip Cover · White Table Base Color · Custom Table Base Graphics	BRANDING KIT 5 - S1175 · Color Edge Print · Custom Table Base Color · Custom Table Base Graphics
BRANDING KIT 3 - S1173 · Deboss Logo · Custom Table Base Color · Custom Table Base Graphics	BRANDING KIT 6 - S1176 · Color Edge Print · White Table Base Color · Custom Table Base Graphics

CUSTOM BASE OPTIONS:

S1178 - Table Base Graphics
PT9020-CC - Pick Custom Base Color
S2000 - Custom Table Base Powder Coat Color and Table Base Graphics



STANDARD BASE COLORS:



White


CUSTOM BASE COLORS:

	Black
	Light Gray
	Charcoal
	Bright Blue
	Royal Blue
	Navy Blue
	Red
	Sky Blue

	Burgundy
	Green
	Forest Green
	Purple
	Yellow
	Orange
	Brown

CUSTOM CUSHION BRANDING:

PT3027 - Logo Slip Cover (27"W x 12"H) PT3030 - Logo Slip Cover (30"W x 12"H) A4451 - Logo Slip Cover (PT1104-SL) (30"W x 12"H)	
LOGO CP - Color Edge Print (12"W x 12"H)	

LOGO DB - Deboss Logo (12"W x 12"H)	
DEBPLT12 - Deboss Plate	

COLOR EDGE UPHOLSTERY CHOICES:

	Blue Ribbon		Celestial
	Purple Passion		Valentine
	Crimson		Dolphin
	Black		

**Color Edge Print upholstery colors do not match Standard Upholstery colors.