

HIP & KNEE FLEXION TREATMENT TABLE

STANDARD FEATURES:

- Load capacity: 800 lbs.
- Solid hardwood construction
- Dual brace support
- 2" High density foam cushion
- Rounded corners
- Manual lift back
- Elevating back: 0 - 70°
- Elevating legs: 0 - 90°
- Hip/Knee flexion
- Bottom shelf for added storage



PT9080



Manual Lift Back



Hip/Knee Flexion

STANDARD TABLES:

PT9080

30"W x 81.5"L x 32"H

STANDARD UPHOLSTERY CHOICES:




STANDARD WOOD CHOICES:



WOOD STAIN OPTIONS:

ESS-Stain

CUSTOM STAIN CHOICES:



CUSTOM CUSHION BRANDING:


TIS 30 COVER/LOGO

Logo Slip Cover
(Lift Back Only)
(30"W x 12"H)



Logo DB


Deboss Logo
(12"W x 12"H)
DEBPLT12 - Deboss Plate



Logo CP - Color Edge Print
(12"W x 12"H)



COLOR EDGE UPHOLSTERY CHOICES:



**Color Edge Print upholstery colors do not match Standard Upholstery colors.

