

WALL MOUNT FOLDING TREATMENT TABLE

STANDARD FEATURES:

- Load capacity: 450 lbs.
- Solid hardwood frame
- 2" High density foam cushion
- Durable wall mounting plate included



WMFTT-2472

STANDARD TABLES:

WMFTT-2472

24"W x 72"L x 31"H

OPTION:

NH - Nose Hole Cushion (on left side only)

STANDARD UPHOLSTERY CHOICES:

 Royal Blue	 American Beauty
 Imperial Blue	 Lark
 Blue Ridge	 Burgundy
 Deep Violet	 Dove
 Forest	 Charcoal
 Mandarin Orange	 Mocha
 Black	 Putty


STANDARD WOOD CHOICES:

 Oak
 Maple

WOOD STAIN OPTIONS:

ESS-Stain

STAIN CHOICES:

 Warm Copper	 Rustic Fall
 Deep Garnet	 Warm Caramel
 Fossil Gray	 Black Beauty
 Tungsten	