

# simpleABI<sup>with CUFF-LINK</sup> 400CL

**Automated Doppler-Based Single-Level ABI + TBI Testing**

**Dual-Testing Capabilities:** Includes Doppler-Free PVR-Based Rapid-Test

Newman Medical

## INDUSTRY-LEADING YEAR WARRANTY

2

ALL PARTS & LABOR  
ACCIDENTAL BREAKAGE  
NO DIAGNOSTIC FEES  
NO HIDDEN FEES  
MINIMAL DOWN TIME

MADE IN  
  
USA

## REAL PEOPLE REAL SUPPORT

NO ROBOTS.  
NO RUNAROUND.  
JUST REAL HELP.



The simpleABI-400CL automated Cuff-Link system is a sophisticated solution for clinics seeking efficient PAD testing. Offering dual testing capabilities, it performs both the Doppler-free, rapid ABI-Q exam (allowing any staff member to quickly screen patients) and comprehensive, Doppler-based ABI and TBI exams. Integrating automated cuff selection and inflation with these advanced capabilities streamlines workflows, minimizes manual steps, and boosts diagnostic accuracy, freeing clinicians to focus on patient care and precise interpretation.

## simpleABI-400CL Features

- Made in the USA with high-quality materials
- Comprehensive 2-year warranty coverage
- Expert support: on-site, online and phone assistance
- Dual testing: ABI-Q rapid-test and Doppler-based exams.
- Automated cuff selection and inflation to reduce user error.
- Versatile testing capabilities for various patient needs.
- Accurate and reliable results.
- User-friendly interface with minimal training required.
- Streamlined workflow

## Ideal for:

- Primary care clinics
- Podiatry offices
- Wound care centers
- Community Health



## simpleABI-400CL Exams



Rapid ABI-Q exam



ABI-Q plus exercise  
(no gym equipment needed)



Doppler ABI



TBI

## simpleABI-400CL Includes:

- simpleABI software
- 4-port Cuff-Link with remote
- 8 MHz Doppler probe
- DPPG probe
- Cuffs
  - 2 x 10cm EZ-Clean Vasocuffs
  - 2 x 12cm EZ-Clean Vasocuffs
  - 2 x 1.9cm disposable digit cuffs
- Adjustable-height roll stand
- User manual / Quick reference guide
- Two-year warranty
  - Includes all parts, labor, and accidental breakage

## simpleABI-400CL Accessories:

### SKU

ABI-120  
STND-500  
ACC-190  
D5  
D8  
SPPG  
SFTW-100  
ABI-WARR

### Description

Microsoft® Windows® PC  
Variable-height roll stand for Cuff-Link systems  
Large simpleABI carry bag  
5MHz vascular probe only  
8 MHz vascular probe only  
Cuff-Link PPG probe - includes toe cuff  
DICOM software integration  
1 year additional warranty for simpleABI systems

## simpleABI-400CL Order Info:

SKU	Description
ABI-400CL	simpleABI-400CL Cuff-Link system