

simpleABI^{with CUFF-LINK} 450CL

Automated Doppler-Based Single-Level ABI + TBI Testing with Stress/Exercise

Dual-Testing Capabilities: Includes Doppler-Free PVR-Based Rapid-Test

INDUSTRY-LEADING YEAR WARRANTY

2

ALL PARTS & LABOR
ACCIDENTAL BREAKAGE
NO DIAGNOSTIC FEES
NO HIDDEN FEES
MINIMAL DOWNTIME



REAL PEOPLE REAL SUPPORT

NO ROBOTS.
NO RUNAROUND.
JUST REAL HELP.



The simpleABI-450CL automated Cuff-Link system offers a versatile and efficient solution for comprehensive PAD assessment, going beyond basic ABI testing. Combining the convenience of automated cuff inflation with a broad range of exam capabilities, including exercise and stress testing, this system excels at uncovering hidden PAD and providing more accurate diagnoses. Its dual testing capability includes the rapid, Doppler-free ABI-Q for quick screening by any staff member, as well as comprehensive Doppler-based exams for in-depth analysis.

simpleABI-450CL Features

- Made in the USA with high-quality materials
- Comprehensive 2-year warranty coverage
- Expert support: on-site, online and phone assistance
- Dual testing: ABI-Q rapid-test and Doppler-based exams.
- Comprehensive PAD testing with exercise/stress testing.
- Automated cuff selection and inflation to reduce user error.
- Accurate diagnoses.
- Improved workflow.

Ideal for:

- Primary care clinics
- Podiatry offices
- Wound care centers
- Community Health



simpleABI-450CL Exams



Rapid ABI-Q exam



ABI-Q plus exercise
(no gym equipment needed)



Doppler ABI



TBI



Stress/Exercise

simpleABI-450CL Includes:

- simpleABI software
- 4-port Cuff-Link with remote
- 8 MHz Doppler probe
- SPPG probe
- Cuffs
 - 2 x 10cm EZ-Clean Vasocuffs
 - 2 x 12cm EZ-Clean Vasocuffs
 - 2 x 1.9cm disposable digit cuffs
- Adjustable-height roll stand
- User manual / Quick reference guide
- Two-year warranty
 - Includes all parts, labor, and accidental breakage

simpleABI-450CL Accessories:

SKU

ABI-120
STND-500
ACC-190
D5
D8
SPPG
SFTW-100
ABI-WARR

Description

Microsoft® Windows® PC
Variable-height roll stand for Cuff-Link systems
Large simpleABI carry bag
5MHz vascular probe only
8 MHz vascular probe only
Cuff-Link PPG probe - includes toe cuff
DICOM software integration
1 year additional warranty for simpleABI systems

simpleABI-450CL Order Info:

| SKU | Description |
|-----------|----------------------------------|
| ABI-450CL | simpleABI-450CL Cuff-Link system |