

# BASIC PRO

## PORTABLE CHIROPRACTIC TABLE

### STANDARD FEATURES:

- Load capacity: 400 lbs.
- Table weight: 46 lbs.
- Height range: 18.5-22.4"
- Folded table: 21.5"W x 37"L x 8"H
- Thoracic & pelvic drops
- Narrow shoulder design
- 5 Angle adjustable headpiece
- Adjustable ankle extension
- Adjustable height, folding legs
- 2" Multi-layer foam cushion
- Black upholstery
- Hand rests
- Paper holder
- Black carry case



BASIC PRO

### STANDARD TABLE:

#### BASIC PRO

27.5"W x 70"L x 18.5-22.4"H



Adjustable Height Folding Legs



Height Range: 18.5-22.4"



Adjustable Headpiece