



AfterMarket
Parts • Service • Support • Accessories

Service Parts List

Click on the model you need parts for

Overhead Stirrer e-A51 Series

[Document update](#)

[e-A51ST020](#)

[e-A51ST040](#)








[e-A51ST060](#)


[e-A51ST100](#)









[e-A51ST200](#)

Note: In one box models are using same spare part

Date	Update
10/20/2020	Initial release
5/6/2021	Replace the incorrect Material P/N listed in Spare part list for e-A51 as per Farrell (OHCT)
12/11/2024	Replace incorrect Material Part number listed in SPL per Lauren S (TRUS)

Order P/N	Material P/N	Description	Image
30586810	40002654	PCB main 230V 60 e-A51	
30586811	40002841	PCB main 230V 40 e-A51	
30586813	40002843	PCB main 230V 20 e-A51	
30586815	40002655	PCB main 115V 60 e-A51	
30586816	40002842	PCB main 115V 40 e-A51	
30586817	40002844	PCB main 115V 20 e-A51	
30586819	40002652	PCB main 230V 100 e-A51	
30586820	40002653	PCB main 115V 100 e-A51	
30586822	40002650	PCB main 230V 200 e-A51	
30586823	40002651	PCB main 115V 200 e-A51	
30586790	10000870	Shaft support e-A51	
30586791	10008309	Knob e-A51	
30586792	10007541	Power Socket e-A51	
30586793	10007566	Power Switch e-A51	
30586794	10007585	Rubber Cover USB e-A51	
30586795	10007586	Rubber Cover Option e-A51	

30586796	10007587	Cylinder shaft e-A51	
30586797	10007589	Rubber Cover shaft e-A51	
30586798	10007590	Rubber Cover screw e-A51	
30586799	10007593	Gasket frontal 200 e-A51	
30586806	10007725	Gasket frontal 100-60-40-20 e-A51	
30586800	10007602	Gasket power socket e-A51	
30586801	10007619	Gasket base e-A51	
30586802	10007648	Rubber Cover cilinder e-A51	
30586803	10007651	Film PCB protection e-A51	

30586804	10007692	Shim driver insulation e-A51	
30586805	10007725	O-ring e-A51	
30586807	40002577	Display e-A51	
30586808	40002606	Chuck assembly e-A51	
30586809	40002731	Collar chuck e-A51	
30586812	40002496	Motor Assembly reduction 60 e-A51	
30586814	40002819	Motor Assembly reduction 20-40 e-A51	
30586818	40002495	Motor Assembly reduction 100 e-A51	
30586821	40002468	Motor Assembly reduction 200 e-A51	
30586824	10008311	Function Label 20 e-A51	
30586825	10008312	Function Label 40 e-A51	
30586826	10008313	Function Label 60 e-A51	
30586827	10008314	Function Label 100 e-A51	
30586828	10008315	Function Label 200 e-A51	
30586836	40002847	Box complete e-A51	
30586847	10007635	Fitting Panel Display e-A51	

30586838	10007180	Power Cable EU e-A51	
30586841	10008269	Power Cable ZH e-A51	
30586768	10007181	Power Cable US e-A51	
30586769	10007397	Power Cable UK e-A51	
30586843	N/A (OHAUS MATERIAL)	Power Cable AU e-A51	

Note: Order P/N is the number for Dealers/ MOs' order; The Material P/N is the raw material number, which might be also printed on the packing label. Here list both part numbers in the part list, in case of any confusion.