

# ESSENTIAL THERA-P™ BARIATRIC ELECTRIC TREATMENT TABLE

## STANDARD FEATURES:


- Load capacity: 650 lbs.
- 34" or 36"W x 77"L x 17-35"H
- Durable and stylish steel base (PET9024/EPT9025)
- White anti-bacterial and anti-microbial powder coat base
- Whisper quiet lift system with single foot control pedal
- Height adjustment: 17-35"
- 34" or 36" Wide cushions in 4 different configuration options
- 2" High density foam cushions with rounded corners
- Individual locking casters
- Head section angle range: +40° to -90°
- Lower extremity angle range: 0 to 65°
- Armrest horizontal range: 0 to 45°
- Armrest vertical range: 0 to 7"
- ADA Compliant



EPT3012

**NOTE:** International (240V) option available

## TABLE CONFIGURATIONS:



**EPT3011**  
1 Section  
34" Width

**EPT3012**  
2 Section  
34" Width

**EPT3013**  
3 Section  
34" Width


**EPT3015**  
5 Section  
34" Width

**EPT3611**  
1 Section  
36" Width

**EPT3612**  
2 Section  
36" Width

**EPT3613**  
3 Section  
36" Width

**EPT3615**  
5 Section  
36" Width




EPT3011 or EPT3611



EPT3012 or EPT3612



EPT3013 or EPT3613



EPT3015 or EPT3615

# OPTIONS:

## OPTIONS:

**EPT9025** - 650 lb Lifting Capacity

**S1136** - Manual Elevating Midsection (only available on EPT3013, EPT3015, EPT3613, & EPT3615)

**S1077** - Dual Foot Control Pedal


**S1078** - Dual Foot Bar Pedal



Manual Elevating Midsection (S1135)



Dual Foot Control Pedal (x2) (S1077)



Dual Foot Bar Control (S1078)

## STANDARD FRAME COLORS:



White

## CUSTOM BASE COLORS:



## STANDARD UPHOLSTERY CHOICES:



## CUSTOM BASE OPTIONS:

**S1178B** - Table Base Graphics

**PT9024-CC** - Pick Custom Base Color

**S2000B** - Custom Table Base Powder Coat Color and Table Base Graphics

