

simpleABI-600CL Automated Cuff-Link System

Newman Medical



Specialists



Cardiology



Vascular Surgery



Interventional Radiology



Rapid ABI-Q



Single-Level



Multi-Level



TBI w/PPG



Exercise

CLINIC PLACEMENT

EXAMS PERFORMED



American-Made Quality



Includes labor, parts, and accidental breakage



Automatic Cuff Selector



Push-Button Remote



simpleABI Roll Stand



ABI-Q Rapid Screen



DICOM Compatible

simpleABI with CUFF-LINK

The simpleABI-600CL Automated Cuff-Link System offers accurate and efficient PAD testing for confident, early diagnosis. With versatile testing options, streamlined workflows, and simplified reimbursement processes, simpleABI enhances convenience and safety for both clinicians and patients. Backed by a 2-year warranty, it's the smart choice for proactive PAD management and improved patient outcomes.

ABI-600CL Exams

- ABI-Q Rapid Screen
- Single-Level
- Multi-Level
- TBI with PPG
- Exercise & Stress

ABI-600CL Includes:

- simpleABI software
- 4-port Cuff-Link with remote
- 8 MHz Doppler probe
- DPPG probe
- Cuffs
 - 2 x 10cm EZ-Clean Vasocuffs
 - 4 x 12cm EZ-Clean Vasocuffs
 - 2 x 12cm long EZ-Clean Vasocuffs
 - 2 x 1.9cm disposable digit cuffs
- Adjustable-height roll stand
- User manual / Quick reference guide
- Two-Year Warranty
 - Includes all parts, labor, and accidental breakage



simpleABI Vascular Systems: Quality You Can Trust

Made in the USA

- High-quality components, manufactured in Colorado.
- Durable design for busy clinical settings.
- 2-year warranty with complete coverage.

Newman Medical: Your Partner in Vascular Care

- Leaders in ultrasound diagnostics since the 1970s.
- Committed to innovation and providing exceptional vascular products.
- Focused on service and long-term customer relationships.

First-Class Support

- Experienced customer service team with extensive industry knowledge.
- Dedicated field representatives for personalized assistance.
- Comprehensive clinical support resources, both online and in-person.



ABI-600CL



Exams Performed



ABI-Q Rapid Screen



Single-Level



Segmental



TBI with PPG



Exercise Exams

Ordering Information

System

SKU

ABI-600CL

DESCRIPTION

simpleABI Multi-Level with Exercise Cuff-Link System

Accessories

SKU

ABI-120
STND-500
ACC-190
D5
SPPG
SFTW-100
ABI-WARR

DESCRIPTION

Microsoft® Windows® PC for use with simpleABI Systems
Variable-Height Roll Stand for Cuff-Link Systems
Large simpleABI carry bag
5Mhz vascular probe only
PPG Probe for use with Cuff-Link Systems - includes toe cuff
DICOM software integration
1 year additional warranty for simpleABI systems

Replacement Cuffs

SKU

CUFF-200
CUFF-210
CUFF-210L
CUFF-230
CUFF-250
CUFF-220D
CUFF-240D
CUFF-SET-242

DESCRIPTION

EZ Clean Vascular Cuff, 10cm (ankle, small arm)
EZ Clean Vascular Cuff, 12cm (arm, ankle, calf)
EZ Clean Vascular Cuff, 12cm long (above knee)
EZ-Clean Vascular Cuff, 17cm Contoured (thigh)
EZ-Clean Vascular Cuff, 7cm long (metatarsal)
Digit Cuff, 2.5cm - Latex-Free
Digit Cuff, 1.9cm - Latex-Free
EZ Clean Vascular Cuff Set - 2x10cm, 4x12cm, 2x12cm - 8/pk



ABI-600CL Technical Specifications

Height: 107 cm (42 inches)

Weight: 5.5 kgs (12 pounds)

Level of Protection against electrical shock: Type B Applied Part, Class II Equipment

Comply with the following Standards: IEC60101-1, IEC60601-2, IEC60601-2-37

Operating Temperature: 10° ~ 40°C (50° ~ 104°F)

Operating Humidity: 30% ~ 75%

Transport/Storage Temperature: -20° ~ 50°C (-4° ~ 122°F)

Transport/Storage Humidity: 5% ~ 90%, non-condensing

Battery Life: 10 hours

GTIN: 852785008609

REIMBURSEMENT INFORMATION

Newman Medical's simpleABI vascular systems frequently utilize CPT codes 93922, 93923, and 93924 for billing and reimbursement purposes.

To obtain expert guidance on the appropriateness of these codes for specific systems, we highly recommend consulting with one of our experienced vascular specialists.

