

# THERA-P™ BOBATH

## ELECTRIC TREATMENT TABLE

### STANDARD FEATURES:


- Load capacity: 650 lbs.
- 40"W x 77"L x 17-35"H
- Durable and stylish steel base (PT9024/PT9025)
- White anti-bacterial and anti-microbial powder coat base
- Whisper quiet lift system with single foot control pedal
- Height adjustment: 17-35"
- 40" Wide cushions in 4 different configuration options
- 2" High density foam cushions with rounded corners
- Low profile locking casters
- Head section angle range: +40° to -90°
- Lower extremity angle range: 0 to 65°
- Armrest horizontal range: 0 to 45°
- Armrest vertical range: 0 to 7"
- ADA Compliant



PT1112

**NOTE:** International (240V) option available

### TABLE CONFIGURATIONS:



1 Section

2 Section

3 Section


5 Section

**PT1111**  
1 Section

**PT1112**  
2 Section

**PT1113**  
3 Section

**PT1115**  
5 Section



PT1111

PT1112

PT1113

PT1115

# OPTIONS:

## OPTIONS:

**EPT9025** - 650 lb Lifting Capacity

**S1136** - Manual Elevating Midsection (only available on PT1113, PT1115)

**S1077** - Dual Foot Control Pedal


**S1078** - Dual Foot Bar Pedal



Manual Elevating Midsection (S1136)



Dual Foot Control Pedal (x2) (S1077)



Dual Foot Bar Control (S1078)

## STANDARD BASE COLORS:



White

## CUSTOM BASE COLORS:



## STANDARD UPHOLSTERY CHOICES:



## CUSTOM BASE OPTIONS:

**S1178B** - Table Base Graphics

**PT9024-CC** - Pick Custom Base Color

**S2000B** - Custom Table Base Powder Coat Color and Table Base Graphics

