

NEW!



# Transferpette® pro

A New Twist on Comfort and Control  
Simple. Easy. Efficient

[BRANDTECH.com](http://BRANDTECH.com)

## Transferpette® pro micropipette

The addition of the Transferpette® pro to the trusted Transferpette® family offers scientists even more choice in selecting the ideal pipette for their unique needs. Designed for repeatability and comfort, the Transferpette® pro is built to meet the demands of your daily lab routine.

The Transferpette® pro features a slim, ergonomic handle paired with a rotating finger rest—a unique design that adjusts to fit any hand position, left- or right-handed. Combined with two intuitive calibration systems—Easy Calibration for tool-free, permanent adjustment and User Adjustment for quick, reversible fine-tuning—the Transferpette® pro ensures precise performance, even with changing liquid types, temperatures, or style of tips.

## Your benefits:

- + **Simple:** Comfortable pipetting with a slim profile and adjustable finger rest
- + **Easy:** Adapt quickly to changing conditions with new user adjustment
- + **Efficient:** Reliable, durable, and convenient

# A New Twist on Comfort and Control



Every application has different requirements. The Transferpette® *pro* adapts perfectly to your individual needs with additional innovative features:

## Adjustable finger rest



Adjustable to 9 positions for an ergonomic hand position

## User Adjustment



Easy and temporary adjustment for special conditions such as temperature, tip shape, or liquid type

## Counter protection



Protects the pipette by preventing over-rotation, ensuring lasting performance and reliability

## Extended volume range



Now available in 300 µl, 2.5 ml, and 1250 µl nominal volume

## Counter protection

Prevents damage from over-rotating beyond volume range

## Rotatable finger rest

Simple adjustment without tools

## User Adjustment

Temporary adjustment to changing conditions

## Extended volume range

13 single-channel pipettes from 0.1 µl to 10 ml

6 multi-channel pipettes from 0.5 µl to 1250 µl



## Transferpette® *pro* and Transferpette® *S* in comparison

	Transferpette® <i>pro</i>	Transferpette® <i>S</i>
Volume range single-channel (adjustable)	12 sizes: from 0.1 µl to 10 ml	10 sizes: from 0.1 µl to 10 ml
Volume range multi-channel	6 sizes: from 0.5 µl to 1250 µl	5 sizes: from 0.5 µl to 300 µl
Size of handle	Slim	Standard
Adjustable finger rest	✓	-
One-handed operation	✓	✓
Uses universal tips	✓	✓
Volume lock	✓	✓
Counter protection	✓	-
Easy Calibration Technology	✓	✓
Temporary User Adjustment	✓	-
Completely autoclavable (121 °C)	✓	✓
UV-resistant	✓	✓
Warranty for single-channel pipettes	3 years	3 years
Warranty for multi-channel pipettes	1 year	1 year

