

# THERA-P PELVIC HEALTH TREATMENT TABLE

## STANDARD FEATURES:

- Load capacity: 650 lbs.
- 30"W x 77"L x 17-35"H
- Durable and stylish steel base with white powder coat finish
- Whisper quiet lift system with single foot control pedal
- Electric elevating midsection
- Single, sliding, adjustable stabilization arm - removable for use on either side of table
- Height adjustment: 17-35"
- 2" High density foam cushions with rounded corners
- Low profile locking casters
- Head section angle range: +40° to -90°
- Lower extremity angle range: 0 to 65°
- Armrest horizontal range: 0 to 45°
- Armrest vertical range: 0 to 7"
- 5 Section configuration includes electric elevating midsection with hand control
- ADA Compliant



PT9045

## STANDARD TABLE:

### PT9045

30"W x 77"L x 17-35"H

QTY: \_\_\_\_\_

## STANDARD UPHOLSTERY CHOICES:

	Royal Blue		American Beauty
	Imperial Blue		Lark
	Blue Ridge		Burgundy
	Deep Violet		Dove
	Forest		Charcoal
	Mandarin Orange		Mocha
	Black		Putty



PT9045 with Dual Footbar Control, Battery, and Custom Base Color

# OPTIONS:

## OPTIONS:

- S1078** - Dual Foot Bar Control
- S2078** - Dual Foot Bar Control with Front Foot Bar Adder
- S1077** - Dual Foot Control Pedal
- PT10001** - Nose Hole Plug
- S1151** - (1) Battery & (1) Wall Charger
- S1152** - (2) Batteries & (1) Wall Charger
- PHO-1040** - Additional Stabilization Arm

## STANDARD FRAME COLORS:



White

## CUSTOM BASE COLORS:



## BRANDING KITS:

### Branding Kit 1 - S1170

- Slip Cover
- Custom Table Base Color
- Custom Table Base Graphics

### Branding Kit 2 - S1171

- Slip Cover
- White Table Base Color
- Custom Table Base Graphics

### Branding Kit 3 - S1173

- Deboss Logo
- Custom Table Base Color
- Custom Table Base Graphics

### Branding Kit 4 - S1174

- Deboss Logo
- White Table Base Color
- Custom Table Base Graphics

### Branding Kit 5 - S1176

- Color Edge Print
- Custom Table Base Color
- Custom Table Base Graphics

### Branding Kit 6 - S1177

- Color Edge Print
- White Table Base Color
- Custom Table Base Graphics

## CUSHION BRANDING:

### PT3030

Logo Slip Cover  
(30"W x 12"H)



**Logo DB** - Deboss Logo  
(12"W x 12"H) **DEBPLT12** -  
Deboss Plate



**Logo CP** - Color Edge Print  
(12"W x 12"H)



## COLOR EDGE UPHOLSTERY CHOICES:



\*\*Color Edge Print upholstery colors do not match Standard Upholstery colors.

## BASE OPTIONS:

- S1178** - Table Base Graphics
- PT9020-CC** - Pick Custom Base Color
- S2000** - Custom Table Base Powder Coat Color and Table Base Graphics

