

Colored MicroTube Caps

New Single Color Caps for 2mL MicroTubes

- 🔗 Screw Caps with Silicone O-Ring
- 🔗 Natural, Assorted, or Single Color Caps
- 🔗 Natural Color Tubes with No Markings
- 🔗 Ideal for Reagents and Kits



By popular demand for researchers, as well as diagnostic / reagent packaging applications, MTC's O-ring screw caps are now available in single color packs of 1000. Order now in blue, green, orange, red, yellow, or natural color.



Colored MicroTube Caps

New Single Color Caps for 2mL MicroTubes

ClearSeal™ Microtubes are made of 100% natural/ colorless polypropylene. They are label ready, as they have no frosted marking spot and no graduations or markings. A silicone O-ring ensures that the screw cap provides a leak-proof seal throughout the temperature range of -90°C to +121°C.

ClearSeal™ tubes are RNase, DNase, pyrogen, and endotoxin free and may be centrifuged at up to 22,000xg.

Dimensions: 13.1×46.8mm.



Cat No.	Description	Qty
Caps for 2.0mL Microtubes		
C3175-B	Screw Cap for 2.0mL MicroTube, Silicone O-Ring, Blue Color, Non-Sterile	1000 / pk
C3175-G	Screw Cap for 2.0mL MicroTube, Silicone O-Ring, Green Color, Non-Sterile	1000 / pk
C3175-O	Screw Cap for 2.0mL MicroTube, Silicone O-Ring, Orange Color, Non-Sterile	1000 / pk
C3175-R	Screw Cap for 2.0mL MicroTube, Silicone O-Ring, Red Color, Non-Sterile	1000 / pk
C3175-Y	Screw Cap for 2.0mL MicroTube, Silicone O-Ring, Yellow Color, Non-Sterile	1000 / pk
C3175-N	Screw Cap for 2.0mL MicroTube, Silicone O-Ring, Natural Color, Non-Sterile	1000 / pk
C3175-AST	Screw Cap for 2.0mL MicroTube, Silicone O-Ring, Non-Sterile, 5 Single Color Bags of 200	1000 / pk
Capless 2.0mL Microtubes		
C3172	Screw Cap Microtube, 2.0mL, Capless, No Graduations or Printing, Self-Standing, Non-Sterile	1000 / cs