

Short Serological Pipettes

Reduced Length 5, 10, and 25mL Pipettes




- ▶ Reduced length - 230mm / 9in.
- ▶ Fits all major pipette controllers
- ▶ Ideal for use in hoods
- ▶ RNase, DNase, and pyrogen-free
- ▶ Positive + negative graduations
- ▶ Non-cytotoxic
- ▶ Barrier plug protects controller
- ▶ Gamma sterilized to SAL 10⁻⁶




P8270-5RL - 5mL | Blue Coded | 99mm Shorter than Standard 5mL



P8270-10RL - 10mL | Orange Coded | 107mm Shorter than Standard 10mL



P8270-25RL - 25mL | Magenta Coded | 85mm Shorter than Standard 25mL



MTC Bio's short serological pipettes feature a reduced length for improved ergonomics.

The design helps to reduce stress or fatigue and prevents accidental contact with hoods, walls, or benchtop equipment.

Cat No.	Description	Qty
P8270-5RL	Serological Pipette, 5mL, Short, 230mm (9 in.), Sterile, Individually Wrapped, 4 packs of 50	200 / cs
P8270-10RL	Serological Pipette, 10mL, Short, 230mm (9 in.), Sterile, Individually Wrapped, 4 packs of 50	200 / cs
P8270-25RL	Serological Pipette, 25mL, Short, 230mm (9 in.), Sterile, Individually Wrapped, 4 packs of 50	200 / cs