

6 LEG CLASSIC WOOD TREATMENT TABLE

STANDARD FEATURES:

- Load capacity: 800 lbs.
- Solid hardwood 6 leg h-brace construction
- 2" High density foam cushion
- Rounded corners
- Triple-bolted leg assembly
- Bottom shelf for added storage



TESS-6L with Cabinet and Drawers

FLAT TOP TABLES:

FESS-6L - Flat Cushion
30"W x 78"L x 31"H

FESS-6LS - Flat, Split Leg Cushion
30"W x 78"L x 31"H

LIFT BACK TABLES:

TESS-6L - Lift Back Cushion
30"W x 78"L x 31"H

TESS-6LS - Lift Back, Split Leg Cushion
30"W x 78"L x 31"H



Shelf, Cabinet & (2) Drawers with Black Wire Pull Handles



Manual Lift Back with 6 Position Elevation



EPTPH01

OPTIONS:

ESS Cabinet 6L - Shelf and Open Cabinet:

- Head End Foot End

ESS Drawers 6L - Shelf, Cabinet & (2) Drawers with Black Wire Pull Handles

- Head End Foot End

NH - Nose Hole in Cushion

TSFO-FHDR - Flush Mount Handles on Drawers

EPTPH01 - Paper Roll Holder (On Foot End with Paper Cutter on Head End)




STANDARD WOOD CHOICES:



WOOD STAIN OPTIONS:

ESS-Stain

STANDARD STAIN CHOICES:



CUSTOM BRANDING:

Decal Custom Small
Small Logo Decal
(6"W x 6"H)



Decal Custom Large
Large Logo Decal
(12"W x 12"H)



STANDARD UPHOLSTERY CHOICES:



CUSTOM CUSHION BRANDING:

TIS 30 COVER/LOGO
Logo Slip Cover
(Lift Back Only)
(30"W x 12"H)



Logo CP - Color Edge Print
(12"W x 12"H)



Logo DB
Deboss Logo
(12"W x 12"H)
DEBPLT12 - Deboss Plate



COLOR EDGE UPHOLSTERY CHOICES:



****Color Edge Print upholstery colors do not match Standard Upholstery colors.**