

CLASSIC

WOOD TREATMENT TABLE

STANDARD FEATURES:

- Load capacity: 550 lbs.
- Solid oak or maple wood frame
- 4 Leg h-brace construction
- 2" High density foam cushion
- Rounded corners
- Lift back has 70° of tilt and six position manual elevation
- Dual-bolted folding leg assembly for convenient shipping



FESS-2772

FLAT TABLES:

FESS-2772 - 27" Flat Cushion
27"W x 72"L x 31"H

FESS-2778 - 27" Flat Cushion
27"W x 78"L x 31"H

FESS-3072 - 30" Flat Cushion
30"W x 72"L x 31"H

FESS-3078 - 30" Flat Cushion
30"W x 78"L x 31"H

LIFT BACK TABLES:

TESS-2772 - 27" Lift Back Cushion
27"W x 72"L x 31"H

TESS-2778 - 27" Lift Back Cushion
27"W x 78"L x 31"H

TESS-3072 - 30" Lift Back Cushion
30"W x 72"L x 31"H

TESS-3078 - 30" Lift Back Cushion
30"W x 78"L x 31"H

OPTIONS:

ESS Bottom Shelf - Bottom Shelf

ESS Cabinet - Shelf and Open Cabinet

Head End Foot End

ESS Drawers - Shelf, Cabinet & (2) Drawers with Black Wire Pull Handles

Head End Foot End

NH - Nose Hole in Cushion

TSFO-FHDR - Flush Mount Handles on Drawers

EPTPH01 - Paper Roll Holder (On Foot End with Paper Cutter on Head End)



EPTPH01

OPTIONS:


STANDARD WOOD CHOICES:



WOOD STAIN OPTIONS:

ESS-Stain

CUSTOM STAIN CHOICES:



TESS-2772 with Shelf, Cabinet, and (2) Drawers

CUSTOM BRANDING:

Decal Custom Small
Small Logo Decal
(6"W x 6"H)



Decal Custom Large
Large Logo Decal
(12"W x 12"H)



STANDARD UPHOLSTERY CHOICES:



CUSTOM CUSHION BRANDING:

TIS 30 COVER/LOGO
Logo Slip Cover
(Lift Back Only)
(30"W x 12"H)




Logo CP - Color Edge Print
(12"W x 12"H)



Logo DB
Deboss Logo
(12"W x 12"H)
DEBPLT12 - Deboss Plate



COLOR EDGE UPHOLSTERY CHOICES:



****Color Edge Print upholstery colors do not match Standard Upholstery colors.**