

NSK WOOD TREATMENT TABLE

STANDARD FEATURES:

- Load capacity: 800 lbs.
- 3/4" Plywood top
- Solid hardwood partial 6 leg h-brace construction
- (2) Shelves on head end for added storage
- 2" High density foam cushion
- Rounded corners
- Triple-bolted leg assembly
- Lift back option has 70° of tilt and six position elevation
- No shin-kick design providing extra room for the practitioner



TESS-NSK100 with Lift Back

STANDARD TABLES:

FESS-NSK100

Flat Cushion
30"W x 78"L x 31"H

TESS-NSK100

Lift Back Cushion
30"W x 78"L x 31"H

FESS-NSK100-S

Flat, Split Leg Cushion
30"W x 78"L x 31"H

FESS-NSK100-S

Lift Back, Split Leg Cushion
30"W x 78"L x 31"H

OPTIONS:

EWT-NH - Nose Hole in Cushion

EPTPH01 - Paper Roll Holder (On Foot End with Paper Cutter on Head End)

STANDARD UPHOLSTERY CHOICES:




STANDARD WOOD CHOICES:



WOOD STAIN OPTIONS:

ESS-Stain

CUSTOM STAIN CHOICES:



CUSTOM BRANDING:

Decal Custom Small

Small Logo Decal
(6"W x 6"H)



Decal Custom Large

Large Logo Decal
(12"W x 12"H)



CUSTOM CUSHION BRANDING:

TIS 30 COVER/LOGO

Logo Slip Cover
(Lift Back Only)
(30"W x 12"H)



Logo CP - Color Edge Print

(12"W x 12"H)




Logo DB

Deboss Logo
(12"W x 12"H)

DEBPLT12 - Deboss Plate



COLOR EDGE UPHOLSTERY CHOICES:



****Color Edge Print upholstery colors do not match Standard Upholstery colors.**