

ARS49ML



83.75" x 55.25" x 31" (H x W x D)

Performance series Pharma-Lab 49 cu.ft. all-refrigerator in stainless steel

Highlights:

Intelligent microprocessor digital temperature controller

Improved temperature accuracy with a separate temperature probe for continuous product readout

Additional sensor reads interior air temperature for superior temperature control

Product Features:

| | |
|---|---|
| Purpose-built design | Designed and engineered to store vaccines and other sensitive materials under stable temperature conditions |
| Adjustable temperature | Interior is pre-set to 5°C, with an adjustable control range from 2 to 10°C |
| Temperature alarm | Visual and audible alarms alert if the interior temperature rises or falls out of set range |
| Complete stainless steel construction | Interior and exterior are entirely constructed from stainless steel for lasting durability and easy cleanup |
| Forced air cooling | Top-mounted commercial fan ensures rapid temperature recovery |
| Automatic defrost with adjustable cycle defrost | Adjustable cycle defrost operation ensures the unit avoid temperature spikes and ice buildup |
| Self-closing doors | Two doors close to save energy but will remain open when the door is opened at a 90° angle for easier loading |
| Open door alarm | Alarm will beep 3 times if the door is left open for 30 seconds, 5 times at 60 seconds, and will sound continuously if the door is open for more than 5 minutes |
| Commercial grade casters | Ships with a set of casters (2 locking) for easy mobility as needed |

ARS49ML Specifications:

| Overview | |
|----------------------------|----------------------|
| Nominal Height | 84.0" (213 cm) |
| Height of Cabinet | 83.75" (213 cm) |
| Nominal Width | 56.0" (142 cm) |
| Width | 55.25" (140 cm) |
| Nominal Depth | 31.0" (79 cm) |
| Depth | 31.0" (79 cm) |
| Depth with Handle | 32.75" (83 cm) |
| Depth with door at 90° | 56.5" (144 cm) |
| Capacity | 49.0 cu.ft. (1388 L) |
| Defrost Type | Automatic |
| Door | Stainless Steel |
| Cabinet | Stainless Steel |
| US Electrical Safety | UL |
| Canadian Electrical Safety | ULC |
| Sanitation | ETL-S |
| Amps | 7.5 |
| Voltage/Frequency | 115 V AC/60 Hz |
| Weight | 410.0 lbs. (186 kg) |
| Parts & Labor Warranty | 2 Years |
| Compressor Warranty | 5 Years |
| UPC | 761101062303 |
| Refrigerator Features | |
| Door Swing | French |
| Shelf Type | Wire |
| Shelf Qty | 8 |
| Adjustable Shelves | Yes |
| Thermostat Type | Digital |
| Fan Type | Interior |
| Refrigerant Type | R290 |
| Refrigerant Amount | 3.2oz. |
| High Side PSI | 290.0 |
| Low Side PSI | 155.0 |
| Wheel Qty | 4 |
| Interior Light | Yes |
| Temperature Range | 2 to 10°C |
| Dimensions | |
| Interior Height | 51.5" (131 cm) |

| | |
|----------------|----------------|
| Interior Width | 50.5" (128 cm) |
| Interior Depth | 27.0" (69 cm) |

