

FF7LW-FS24LSTACKMED2



56" x 24" x 23.5" (H x W x D)

Stacked combination of FF7LWBIMED2 auto defrost all-refrigerator and FS24LMED2 compact manual defrost all-freezer, both with locks, digital controls, and NIST calibrated alarm/thermometers

Highlights:

Specially designed and featured to meet the CDC, VFC and VFA's guidelines for safe vaccine storage

Independent refrigerator and freezer storage in one slim-fitting footprint for space-challenged facilities

High temperature alarm includes NIST calibrated external temperature readout (NIST certificate included)

Product Features:

| | |
|-------------------------------------|---|
| Meets CDC & VFC guidelines | Designed to comply with the most current guidelines for vaccine storage |
| Space-saving design | Allows you to create a full refrigerator-freezer with independent controls in a slim-fitting footprint |
| FF7LWBIMED2 all-refrigerator | 5.5 cu.ft. auto defrost all-refrigerator with digital controls, internal fan, flat door liner, and adjustable shelving |
| FS24LMED2 all-freezer | Compact manual defrost all-freezer with digital thermostat and wire shelving |
| Easy assembly | Unit ships in three boxes for customer to install the StackRack on the bottom refrigerator/freezer and place and secure the other unit on top |
| Digital thermostats | Digital microprocessor for easy and precise temperature management, with all controls externally located for added convenience |
| Access ports | Rear-mounted 3/8" access port in both units allows users to install their own monitoring equipment |
| High/low temperature alarms | Audible alarm sounds if the interior temperature rises or falls out of set range |
| NIST calibrated temperature display | Alarms include continuous NIST calibrated display of the current interior and high/low temperature, and can be switched from Celsius to Fahrenheit at the touch of a button |

FF7LW-FS24LSTACKMED2 Specifications:

| Overview | |
|----------------------------|--------------------|
| Nominal Height | 56.0" (142 cm) |
| Height of Cabinet | 56.0" (142 cm) |
| Nominal Width | 24.0" (61 cm) |
| Width | 24.0" (61 cm) |
| Nominal Depth | 24.0" (61 cm) |
| Depth | 23.5" (60 cm) |
| Capacity | 6.9 cu.ft. (195 L) |
| Door | White |
| Cabinet | White |
| US Electrical Safety | ETL |
| Canadian Electrical Safety | ETL-C |
| Amps | 1.3 |
| Voltage/Frequency | 115 V AC/60 Hz |
| Parts & Labor Warranty | 2 Years |
| Compressor Warranty | 5 Years |
| UPC | 761101073798 |

| Refrigerator-Freezer | |
|-----------------------------------|--------------------|
| Door Swing | RHD |
| Reversible | Factory Reversible |
| Adjustable Shelves | Yes |
| Interior Light | Yes |
| Refrigerator - Interior Height | 26.25" (67 cm) |
| Refrigerator - Interior Width | 21.0" (53 cm) |
| Refrigerator - Interior Depth | 18.0" (46 cm) |
| Refrigerator - Capacity | 5.5 cu.ft |
| Refrigerator - Shelf Type | Wire |
| Refrigerator - Shelf Qty | 3 |
| Refrigerator - Temperature Range | 2° to 8°C |
| Refrigerator - Refrigerant Type | R600a |
| Refrigerator - Refrigerant Amount | 0.64 oz. |
| Freezer - Interior Height | 14.38" (37 cm) |
| Freezer - Interior Width | 13.5" (34 cm) |
| Freezer - Interior Depth | 13.25" (34 cm) |
| Freezer - Capacity | 1.4 cu.ft |

| | |
|------------------------------|-----------------------|
| Freezer - Shelf Type | Wire |
| Freezer - Shelf Qty | 1 |
| Thermostat Type | Digital |
| Fan Type | Interior and Exterior |
| Refrigerant Type | R600a |
| Refrigerant Amount | 0.64 oz. |
| High Side PSI | 180.0 |
| Low Side PSI | 80.0 |
| Freezer - Temperature Range | -20° to -15°C |
| Freezer - Refrigerant Type | R600a |
| Freezer - Refrigerant Amount | 0.8 oz. |

