

FF7LWBIVAC



33.75" x 23.63" x 23.75" (H x W x D)

Built-in undercounter medical all-refrigerator for temperature stable medical storage, with interior basket drawers, internal fan, lock, and more

Highlights:

Medical grade design conforms with CDC guidelines for vaccine storage

Interior drawers ensure stored items are not exposed to room temperature when the door is open

High/low temperature alarm includes NIST calibrated temperature readout of the interior

Product Features:

Improved temperature stability	Interior is designed to operate within 2°C of set temperature performance
Interior basket drawers	Three slide-out wire baskets with plastic fronts offer convenient storage solutions while protecting stored items when the door is open
Digital thermostat	Electronic controls ensure precise and intuitive temperature management
High/low temperature alarm	Audible alarm sounds if the interior temperature rises or falls out of set range
NIST calibrated temperature display	Alarm includes NIST calibrated display of the current interior and high/low temperature
Self-closing door	Door closes automatically when opened at angles less than 90°
3/8" access port	Unit comes with a factory installed probe hole for users to connect their own monitoring equipment
Hospital grade cord and plug	Secure three-pronged 'green dot' cord for added unit safety in high traffic settings
Automatic defrost	Reduced user maintenance with automatic defrost

FF7LWBIVAC Specifications:

Overview	
Nominal Height	34.0" (86 cm)
Height of Cabinet	33.75" (86 cm)
Nominal Width	24.0" (61 cm)
Width	23.63" (60 cm)
Nominal Depth	24.0" (61 cm)
Depth	23.75" (60 cm)
Depth with Handle	25.63" (65 cm)
Depth with door at 90°	46.25" (117 cm)
Capacity	2.6 cu.ft. (74 L)
Defrost Type	Automatic
Door	White
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	1.4
Voltage/Frequency	115 V AC/60 Hz
Weight	100.0 lbs. (45 kg)
Parts & Labor Warranty	2 Years
Compressor Warranty	5 Years
UPC	761101072555
Refrigerator Features	
Door Swing	RHD
Reversible	Factory Reversible
Height Adjustable	1
Thermostat Type	Digital
Fan Type	Interior and Exterior
Refrigerant Type	R600a
Refrigerant Amount	0.64oz.
High Side PSI	180.0
Low Side PSI	80.0
Level Legs Qty	4
Interior Drawers	Yes
Interior Drawers Qty	3
Temperature Range	2° to 8°C
Dimensions	
Interior Height	22.5" (57 cm)
Interior Width	17.25" (44 cm)

Interior Depth 15.0"

