

FF7LWBIVACADA



32.25" x 23.63" x 23.25" (H x W x D)

Built-in ADA undercounter medical all-refrigerator for temperature stable medical storage, with interior basket drawers, internal fan, lock, and more

Highlights:

Medical grade all-refrigerator ideally suited for vaccine and pharmaceutical storage

Interior drawers ensure stored items are not exposed to room temperature when the door is open

Automatic defrost with an internal fan for stabilized cooling designed to stay within 2°C of set temperature

Product Features:

Improved temperature stability	Interior is designed to operate within 2°C of set temperature performance
ADA Compliant	32" height is designed to fit under lower ADA compliant counters
Interior basket drawers	Three slide-out wire baskets with plastic fronts offer convenient storage solutions while protecting stored items when the door is open
Digital thermostat	Electronic controls ensure precise and intuitive temperature management
High/low temperature alarm	Audible alarm sounds if the interior temperature rises or falls out of set range
NIST calibrated temperature display	Alarm includes NIST calibrated display of the current interior and high/low temperature
Self-closing door	Door closes automatically when opened at angles less than 90°
3/8" access port	Unit comes with a factory installed probe hole for users to connect their own monitoring equipment
Hospital grade cord and plug	Secure three-pronged 'green dot' cord for added unit safety in high traffic settings

